[Gebrauchsanleitung](#) | [Operating manual](#) | [Mode d'emploi](#) |  
[Instrucciones de manejo](#) | [Istruzione](#) | [Instruções de utilização](#) | [操作手册](#) |  
[Руководство по эксплуатации](#) | [사용지침](#) | [Bruksanvisning](#) | [使用説明書](#) |  
[Használati utasítás](#) | [Návod k použití](#) | [Gebruiksaanwijzing](#) |  
[Instrukcja użytkowania](#) | [Kullanım Talimatları](#)

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# Table of Contents

<b>1 Introduction .....</b>	<b>49</b>	<b>16 Calibration service .....</b>	<b>90</b>
1.1 Scope of delivery .....	49	<b>17 Warranty .....</b>	<b>91</b>
1.2 Terms of use .....	49	<b>18 Disposal .....</b>	<b>92</b>
<b>2 Safety Instructions .....</b>	<b>50</b>		
2.1 General safety instructions.....	50		
2.2 Function.....	50		
2.3 Limitations of Use .....	52		
2.4 Operating Limitations.....	52		
2.5 Operating Exclusions .....	52		
2.6 Storage Conditions .....	53		
2.7 Recommended Application Range .	53		
2.8 Dispenser Selection Chart .....	54		
<b>3 Functional and operating elements .....</b>	<b>56</b>		
<b>4 Assembly .....</b>	<b>59</b>		
4.1 First Steps .....	59		
4.2 Priming .....	60		
<b>5 Operation .....</b>	<b>62</b>		
5.1 Dispensing.....	62		
5.2 Accessories .....	63		
<b>6 Error limits .....</b>	<b>66</b>		
<b>7 Checking the Volume (Calibration) .....</b>	<b>67</b>		
<b>8 Adjustment.....</b>	<b>68</b>		
8.1 Digital model .....	68		
8.2 Analog-adjustable model .....	69		
<b>9 Cleaning .....</b>	<b>70</b>		
9.1 Cleaning, Analog-adjustable model, Fixed-volume model .....	70		
9.2 Cleaning Digital model .....	71		
9.3 Replacing discharge tube/valves ....	73		
<b>10 Autoclaving .....</b>	<b>75</b>		
<b>11 Troubleshooting .....</b>	<b>76</b>		
<b>12 Product markings .....</b>	<b>78</b>		
<b>13 Ordering Information .....</b>	<b>79</b>		
<b>14 Accessories/spare parts .....</b>	<b>82</b>		
<b>15 Repairs.....</b>	<b>88</b>		
15.1 Sending for repair .....	88		

# 1 Introduction

## 1.1 Scope of delivery

Dispensette® S or Dispensette® S Organic bottle-top dispenser, for GL 45 threaded bottles, dosing cannula or dosing cannula with back-dosing valve, telescopic aspiration tube, back-dosing tube (optional for devices with back-dosing valve), assembly wrench, various bottle adapters, a quality certificate and these instructions for use.

Nominal volume ml	Adapter for bottle thread, PP	Filling tube length mm
1, 2, 5, 10	GL 24-25, GL 28, S 28, GL 32-33, GL 38, S 40	125- 240
25 / 50, 100	GL 32-33, GL 38, S 40	170- 330

## 1.2 Terms of use

- Carefully read the operating manual before using the device for the first time.
- The operating manual is part of the device and must be kept in an easily accessible place.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find up-to-date versions of the operating manual on our website: [www.brand.de](http://www.brand.de).

### 1.2.1 Hazard levels

The following signal words identify possible hazards:

Signal word	Meaning
DANGER	Will lead to serious injury or death.
WARNING	May lead to serious injury or death.
CAUTION	May lead to minor or moderate injuries.
NOTICE	May lead to property damage.

### 1.2.2 Symbols

Symbol	Meaning
	Hazardous area

### 1.2.3 Format

Format	Meaning	Format	Meaning
1. Task	Indicates a task.	>	Indicates a condition.
a., b., c.	Indicates the individual steps of a task.	⇒	Indicates a result.

## 2 Safety Instructions

### 2.1 General safety instructions

#### Please read carefully!

The instrument Dispensette® S can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

1. Every user must read and observe this operating manual before using the device.
2. Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves.
3. Observe all specifications provided by reagent manufacturers.
4. When dispensing inflammable media, make sure to avoid to buildup of static charge, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.
5. Use the instrument only for dispensing liquids, with strict regard to the defined limitations of use and operating limitations. Comply with the operating exclusions (Operating Exclusions, p. 52)! If in doubt, contact the manufacturer or supplier.
6. Always perform work in a manner that does not endanger yourself or other people. When dispensing, the discharge tube must always point away from you or any other person. Avoid splattering. Only use suitable vessels.
7. Never press down the piston when the discharge tube closure is attached.
8. Never remove the discharge tube while the dispensing cylinder is filled.
9. Reagents can accumulate in the screw cap of the discharge tube. Thus, the screw cap should be cleaned regularly.
10. For small bottles, and when using the flexible discharge tube, use a bottle stand to prevent tipping over.
11. Never carry the mounted instrument by the cylinder sleeve or the valve block. Breakage or loosening of the cylinder can lead to personal injury from chemicals (First Steps, p. 59, Fig. 3).
12. Never use force. Use smooth gentle movements to operate the piston upwards and downwards.
13. Use only original accessories and original replacement parts. Do not make any technical modifications. Do not dismantle the instrument any further than is described in the operating manual!
14. Always check that the instrument is in proper working condition before use. The user can come into contact with media if the instrument has been insufficiently cleaned or inspected. If there is a sign of a potential malfunction (e.g., piston difficult to move, sticking valves or leakage), immediately stop dispensing and consult the Troubleshooting section of this manual (Troubleshooting, p. 76). Contact the manufacturer, if necessary.

### 2.2 Function

With the Dispensette® S and Dispensette® S Organic bottle-top dispensers, liquids can be dispensed directly from the supply bottle. Available in Digital, Analog and Fixed models. The instruments are marked DE-M and optionally equipped with recirculation valve.

## 2.2.1 Dispensette S (red color code)

Digital

Analog

Fixed

## 2.2.2 Dispensette S Organic (yellow color code)

Digital

Analog

Fixed

## 2.2.3 Operation

### NOTICE

#### Dispensing hydrofluoric acid (HF)

For dispensing hydrofluoric acid (HF), we recommend the use of the Dispensette® S Trace Analysis bottle-top dispenser with platinum-iridium valve spring (see separate operating manual).

When the instrument is correctly used, the dispensed liquid comes into contact with only the following chemically resistant materials:

#### Dispensette® S

Borosilicate glass, Al<sub>2</sub>O<sub>3</sub>-ceramic, ETFE, FEP, PFA, PTFE, platinum-iridium, PP (screw cap).

#### Dispensette® S Organic

Borosilicate glass, Al<sub>2</sub>O<sub>3</sub>ceramic, ETFE, FEP, PFA, PTFE, tantalum, PP (screw cap).

Alternatively ETFE/PTFE bottle adapters can be used (Accessories/spare parts, p. 82). The suitability of ETFE/PTFE bottle adapters must be checked by the user.

## 2.3 Limitations of Use

This instrument is designed for dispensing liquids, observing the following physical limits:

- Operating temperature from +15 °C to +40 °C (from 59 °F to 104 °F) of instrument and reagent
- Vapor pressure up to max. 600 mbar. Aspirate slowly above 300 mbar, in order to prevent the liquid from boiling
- Kinematic viscosity up to 500 mm<sup>2</sup>/s (dynamic viscosity [mPas] = kinematic viscosity [mm<sup>2</sup>/s] x density [g/cm<sup>3</sup>])
- Density up to 2.2 g/cm<sup>3</sup>

## 2.4 Operating Limitations

- Liquids that form deposits (e.g., crystallizing solutions or concentrated alkaline solutions) may make the piston difficult to move or may cause jamming. If the piston movement becomes sluggish or stiff, the instrument should be cleaned immediately. See also Cleaning, p. 70.
- When dispensing inflammable media, make sure to avoid to buildup of static charge (e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth).
- The instrument is designed for general laboratory applications and complies with the requirements of the relevant standards (e.g., DIN EN ISO 8655). Compatibility of the instrument for a specific application (e.g., trace material analysis or food sector) must be checked by the user. Approvals for specific applications (e.g., for the production and administration of food, pharmaceuticals, or cosmetics) have not been given.

## 2.5 Operating Exclusions

### 2.5.1 Dispensette S

Never use Dispensette® S for:

- Liquids that attack Al<sub>2</sub>O<sub>3</sub>-ceramic, ETFE, FEP, PFA, and PTFE (e.g., dissolved sodium azide\*)
- Liquids that attack borosilicate glass (e.g., hydrofluoric acid)
- Liquids that are decomposed catalytically by platinum-iridium (e.g., H<sub>2</sub>O<sub>2</sub>)
- Hydrochloric acid > 20% and nitric acid > 30%
- Tetrahydrofuran
- Trifluoroacetic acid
- Explosive liquids
- Carbon disulfide
- Suspensions (e.g., activated charcoal) because solid particles may clog or damage the instrument
- Liquids that tend to polymerize spontaneously (e.g., adhesives and their components)
- Liquids that attack PP (screw cap and bottle adapter)\*\*.

\* Dissolved sodium azide permitted up to a concentration of max. 0.1 %.

\*\* Alternatively, ETFE/PTFE bottle adapters and ETFE/PTFE screw caps can be used (accessories from Accessories/spare parts, p. 82). The suitability of ETFE/PTFE bottle adapters and ETFE/PTFE screw caps must be checked by the user.

### 2.5.2 Dispensette S Organic

Never use Dispensette® S Organic for:

- Liquids that attack  $\text{Al}_2\text{-O}_3$ -ceramic, tantalum, ETFE, FEP, PFA, and PTFE (e.g., dissolved sodium azide\*)
- Liquids that attack borosilicate glass (e.g., hydrofluoric acid)
- Bases and salt solutions
- Explosive liquids
- Carbon disulfide
- Suspensions (e.g., activated charcoal) because solid particles may clog or damage the instrument
- Liquids that tend to polymerize spontaneously (e.g., adhesives and their components)
- Liquids that attack PP (screw cap and bottle adapter)\*\*.

\* Dissolved sodium azide permitted up to a concentration of max. 0.1 %.

\*\* Alternatively, ETFE/PTFE bottle adapters and ETFE/PTFE screw caps can be used (accessories from Accessories/spare parts, p. 82). The suitability of ETFE/PTFE bottle adapters and ETFE/PTFE screw caps must be checked by the user.

## 2.6 Storage Conditions

Store the unit and accessories in a cool and dry place in cleaned condition only.

Storage temperature from  $-20$  to  $50^\circ\text{C}$  ( $-4$  to  $122^\circ\text{F}$ ).

## 2.7 Recommended Application Range

Dispensette® S offers a wide range of applications for dispensing aggressive reagents (e.g., concentrated acids such as  $\text{H}_3\text{PO}_4$  and  $\text{H}_2\text{SO}_4$ , bases such as NaOH and KOH, and salt solutions as well as various polar solvents).

Dispensette® S Organic is used for dispensing organic solvents (e.g., chlorinated and fluorinated hydrocarbons such as trichlorotrifluoroethane and dichloromethane or acids such as concentrated HCl and  $\text{HNO}_3$  as well as for trifluoroacetic acid (TFA), tetrahydrofuran (THF) and peroxides).

### NOTICE

#### Dispenser selection

For guidelines on selecting the right dispenser, please observe the corresponding exclusions of use and the following 'Dispenser Selection Chart'.

### NOTICE

#### Dispensing hydrofluoric acid (HF)

For dispensing hydrofluoric acid (HF), we recommend the use of the Dispensette® S Trace Analysis bottle-top dispenser with platinum-iridium valve spring (see separate operating manual).

## 2.8 Dispenser Selection Chart

### NOTICE

#### Hydrofluoric acid (HF)

For dispensing hydrofluoric acid (HF), we recommend the Dispensette S Trace Analysis bottle-top dispenser with platinum-iridium valve spring.

See separate instructions for use at

[www.brand.de/om](http://www.brand.de/om).

	Dispensette® S	Dispensette® S Organic
1,2-diethylbenzene	✓	✓
1,4-dioxane	—	✓
1-butanol	✓	✓
1-decanol	✓	✓
Acetaldehyde	✓	✓
Acetic acid, ≤ 96%	✓	✓
Acetic acid, 100% (glacial acetic acid)	✓	✓
Acetic anhydride	—	✓
Acetone	✓	✓
Acetonitrile	✓	✓
Acetophenone	—	✓
Acetyl chloride	—	✓
Acetylacetone	✓	✓
Adipic acid	✓	—
Allyl alcohol	✓	✓
Aluminum chloride	✓	—
Amino acids	✓	—
Ammonia, ≤ 20%	✓	✓
Ammonia, 20–30%	—	✓
Ammonium chloride	✓	—
Ammonium fluoride	✓	—
Ammonium sulfate	✓	—
Amyl alcohol (pentanol)	✓	✓
Amyl chloride (chloropentane)	—	✓
Aniline	✓	✓
Barium chloride	✓	—
Benzaldehyde	✓	✓
Benzene	✓	✓
Benzoic acid methyl ester	✓	✓
Benzoyl chloride	✓	✓
Benzyl alcohol	✓	✓
Benzyl chloride	✓	✓
Benzylamine	✓	✓
Boric acid, ≤ 10%	✓	✓
Bromobenzene	✓	✓
Bromonaphthalene	✓	✓
Butanediol	✓	✓
Butyl acetate	✓	✓
Butyl methyl ether	✓	✓
Butylamine	✓	✓
Butyric acid	✓	✓

	Dispensette® S	Dispensette® S Organic
Calcium carbonate	✓	—
Calcium chloride	✓	—
Calcium hydroxide	✓	—
Calcium hypochlorite	✓	—
Carbon tetrachloride	—	✓
Chloroacetaldehyde, ≤ 45%	✓	✓
Chloroacetic acid	✓	✓
Chloroacetone	✓	✓
Chlorobenzene	✓	✓
Chlorobutane	✓	✓
Chloroform	—	✓
Chloronaphthalene	✓	✓
Chlorosulfonic acid	—	✓
Chromic acid, ≤ 50%	✓	✓
Chromosulfuric acid	✓	—
Copper sulfate	✓	—
Cresol	—	✓
Cumene (isopropylbenzene)	✓	✓
Cyclohexane	—	✓
Cyclohexanone	✓	✓
Cyclopentane	—	✓
Decane	✓	✓
Dibenzyl ether	✓	✓
Dichloroacetic acid	—	✓
Dichlorobenzene	✓	✓
Dichloroethane	—	✓
Dichloroethylene	—	✓
Dichloromethane	—	✓
Diesel oil (fuel oil), boiling point 250–350°C	—	✓
Diethanolamine	✓	✓
Diethyl ether	—	✓
Diethylamine	✓	✓
Diethylene glycol	✓	✓
Dimethyl sulfoxide (DMSO)	✓	✓
Dimethylaniline	✓	—
Dimethylformamide (DMF)	✓	✓
Diphenyl ether	✓	✓
Essential oils (aromatic oils)	—	✓
Ethanol	✓	✓
Ethanolamine	✓	✓
Ethyl acetate	✓	✓
Ethyl methyl ketone	✓	✓
Ethylbenzene	—	✓
Ethylene chloride	—	✓
Fluoroacetic acid	—	✓
Formaldehyde, ≤ 40%	✓	—
Formamide	✓	✓
Formic acid, ≤ 100%	—	✓
Fuel oil (diesel oil), boiling point 250–350°C	—	✓
Gasoline (petroleum spirit), boiling point 70–180°C	—	✓
Glycerin	✓	✓
Glycol (ethylene glycol)	✓	✓
Glycolic acid, ≤ 50%	✓	—
Heptane	—	✓
Hexane	—	✓

	Dispensette® S	Dispensette® S Organic
Hexanoic acid	✓	✓
Hexanol	✓	✓
Hydrobromic acid	—	✓
Hydrochloric acid, ≤ 20%	✓	✓
Hydrochloric acid, 20–37% **	—	✓
Hydrogen peroxide, ≤ 35%	—	✓
Hydroiodic acid, ≤ 57% **	✓	✓
Isoamyl alcohol	✓	✓
Isobutanol	✓	✓
Isooctane	—	✓
Isopropanol (2-propanol)	✓	✓
Isopropyl ether	✓	✓
Lactic acid	✓	—
Methanol	✓	✓
Methoxybenzene	✓	✓
Methyl butyl ether	✓	✓
Methyl formate	✓	✓
Methyl propyl ketone	✓	✓
Methylene chloride	—	✓
Mineral oil (engine oil)	✓	✓
Monochloroacetic acid	✓	✓
n-amyl acetate	✓	✓
Nitric acid, ≤ 30%	✓	✓
Nitric acid, 30–70% */**/**	—	✓
Nitrobenzene	✓	✓
Oleic acid	✓	✓
Oxalic acid	✓	—
Pentane	—	✓
Peracetic acid	—	✓
Perchloric acid	✓	✓
Perchloroethylene	—	✓
Petroleum ether, boiling point 40–70°C	—	✓
Petroleum, boiling point 180–220°C	—	✓
Phenol	✓	✓
Phenylethanol	✓	✓
Phenylhydrazine	✓	✓
Phosphoric acid, 85	✓	✓
Piperidine	✓	✓
Potassium chloride	✓	—
Potassium dichromate	✓	—
Potassium hydroxide	✓	—
Potassium permanganate	✓	—
Propionic acid	✓	✓
Propylene glycol (propanediol)	✓	✓
Pyridine	✓	✓
Pyruvic acid	✓	✓
Salicylaldehyde	✓	✓
Scintillation cocktail	✓	✓
Silver acetate	✓	—
Silver nitrate	✓	—

	Dispensette® S	Dispensette® S Organic
Sodium acetate	✓	—
Sodium chloride	✓	—
Sodium dichromate	✓	—
Sodium fluoride	✓	—
Sodium hydroxide, ≤ 30%	✓	—
Sodium hypochlorite	✓	—
Sulfuric acid, ≤ 98%	✓	✓
Tartaric acid	✓	—
Tetrachloroethylene	—	✓
Tetrahydrofuran (THF) */**	—	✓
Tetramethylammonium hydroxide	✓	—
Toluene	—	✓
Trichloroacetic acid	—	✓
Trichlorobenzene	—	✓
Trichloroethane	—	✓
Trichloroethylene	—	✓
Trichlorotrifluoroethane	—	✓
Triethanolamine	✓	✓
Triethylene glycol	✓	✓
Trifluoroacetic acid (TFA)	—	✓
Trifluoroethane	—	✓
Turpentine	—	✓
Urea	✓	—
Xylene	—	✓
Zinc chloride, ≤ 10%	✓	—
Zinc sulfate, ≤ 10%	✓	—

#### Legend:

✓ = The device is suitable for the medium

— = The device is not suitable for the medium

This table has been carefully tested and is based on the most current information available. Always observe the operating manual of the instrument and the specifications provided by the reagent manufacturer. In addition to the chemicals listed above, a number of organic or inorganic saline solutions (e.g. biological buffers), biological detergents as well as media for cell cultures can be dispensed. If you need chemical hazard statements that are not on the list, you are welcome to contact BRAND.

Last updated: 1123/15

\* Use ETFE/PTFE bottle adapter

\*\* Use PTFE sealing ring for valve block

\*\*\* Dispensette® S Organic, Analog-adjustable\*\*\*.

# 3 Functional and operating elements


## Dispensette® S Digital

- 1 Volume adjustment
- 2 Cover plate
- 3 Housing
- 4 Protective sleeve/dispensing cylinder
- 5 Piston unit
- 6 Valve block
- 7 Valve block adapter (GL 45 bottle thread)
- 8 Recirculation tube
- 9 Telescopic filling tube
- 10 Latch
- 11 Safety lock
- 12 Piston mounting nut
- 13 Knob
- 14 Screw cap
- 15 Discharge tube
- 16 Valve lever, recirculation valve

## Dispensette® S Analog

- 1 Cover plate
- 2 Piston seat
- 3 Housing
- 4 Pointer
- 5 Volume adjustment
- 6 Protective sleeve/dispensing cylinder
- 7 Valve block
- 8 Valve block adapter (GL 45 bottle thread)
- 9 Recirculation tube
- 10 Telescopic filling tube
- 11 Screw cap
- 12 Discharge tube
- 13 Valve lever, recirculation valve
- 14 Piston unit

### Telescoping filling tube and recirculation tube

## Assembly tool

## 4 Assembly

### 4.1 First Steps


#### ⚠ WARNING

#### Follow the safety instructions

- Wear protective clothing, eye protection and protective gloves!
- Always wear protective gloves when touching the instrument or the bottle, especially when using dangerous liquids.
- Follow all safety instructions and comply with the limitations of use, see Limitations of Use, p. 52.
- Comply with the operating limitations, see Operating limitations, p. 52.

#### 1. Installing the filling tube/recirculation tube

- a. Adjust the length of the telescoping filling tube to the bottle height and attach it.
  - b. Center and attach the filling tube (part with smaller diameter) carefully to avoid damaging the nozzle.
- ⇒ If a discharge tube with a recirculation valve is used, the recirculation tube must also be installed.
- c. Insert the recirculation tube with the opening pointing outward.

#### 2. Mounting and aligning the instrument on a bottle

#### NOTICE

To avoid tipping over, use a bottle stand for small bottles.

- a. Screw the instrument (GL 45 threads) onto the reagent bottle, and then align the discharge tube with the bottle label. This is done by rotating the valve block with the discharge tube.

### 3. Transporting the instrument

#### NOTICE

The scope of delivery includes bottle adapters made of polypropylene (PP). These may be used only for media that do not attack PP. Alternatively, ETFE/PTFE bottle adapters can be used (Accessories/ spare parts, p. 82). The user must review carefully to ensure the suitability of ETFE/PTFE bottle adapters.

- a. For bottles with other thread sizes, select a suitable bottle adapter.
- b. Always wear protective gloves when touching the instrument or the bottle, especially when using dangerous liquids.
- c. When mounted to a reagent bottle, always carry the instrument as shown in the figure.

#### NOTICE

##### Template for adapter

An adapter template is available for download at [www.brand.de](http://www.brand.de) in the Service section. This can be used to determine thread and adapter sizes. This adapter template for the respective product can also be found in the on-line shop <https://shop.brand.de> in the “Selection Charts” section.

## 4.2 Priming

#### ⚠ WARNING

##### Follow during every use, especially with hazardous media

- > Wear protective clothing, eye protection and appropriate hand protection!
- > Never press down the piston when the screw cap is screwed on!
- > Avoid splashing the reagent!
- > To avoid splashes dispense slowly.
- > Liquid may accumulate in the screw cap. Open the screw cap slowly in order to prevent splashing.
- > Follow all safety instructions and comply with the operating exclusions and limitations of use; see Limitations of use, p. 52 and Operating exclusions, p. 52.

#### NOTICE

Before using the instrument for the first time, ensure it is rinsed carefully and discard the first few samples dispensed. Avoid splashes, as ethanol and glycerin residues may be present in the instrument.

## Instruments with recirculation valve

- a. Open the screw cap of the dispensing tube.

- b. Turn valve to 'Recirculate'.

- c. For deaerating gently pull up the piston approx. 30 mm and push it down rapidly until the lower stop. Repeat this process at least 5 times.

- d. Turn valve to 'Dispense'.

- e. To avoid splashes when deaerating hold the discharge tube on the inner wall of a suitable receiving vessel and dispense liquid to vent deaerate the discharge tube until it is bubble-free. Wipe away any remaining drops from the discharge tube.

## Instruments without recirculation valve

- Open the screw cap of the discharge tube (see 'instrument with recirculation valve', Fig. a.). To avoid splashes, hold discharge tube orifice on the inner wall of a suitable receiving vessel.
- For deaerating pull up the piston approx. 30 mm and push it down rapidly until the lower stop. Repeat this procedure approximately 5 times until the discharge tube is bubble-free.

# 5 Operation

## 5.1 Dispensing

### 1. Selecting the volume

**Digital:** Rotate the volume-setting wheel until the desired volume is indicated (mechanical counter).

**Analog:** Loosen the volume selector thumb screw  $\frac{3}{4}$  turn (1), set the pointer to the desired volume (2) and then retighten the volume thumb screw (3).

**Fixed:** The volume is non-adjustable and cannot be changed.

### 2. Dispensing

#### ⚠ WARNING

**Follow during every use, especially with hazardous media**

- > Wear protective clothing, eye protection and appropriate hand protection!
- > Never press down the piston when the screw cap is screwed on!
- > Avoid splashing the reagent!
- > To avoid splashes dispense slowly.
- > Liquid may accumulate in the screw cap. Open the screw cap slowly in order to prevent splashing.
- > Follow all safety instructions and comply with the operating exclusions and limitations of use; see Limitations of use, p. 52 and Operating exclusions, p. 52.

- a. Unscrew the screw cap of the dispensing tube.
- b. When using instruments equipped with a recirculation valve, turn the valve to 'Dispensing'.
- c. Hold the discharge tube orifice on the inner wall of a suitable receiving vessel.

- d. Gently lift the piston until the upper stop and then depress piston slowly and steadily with minimal force until the lower stop.

- e. Wipe off the discharge tube against the inner wall of the receiving vessel.  
f. Seal the dispensing tube with the screw cap.

### NOTICE

After using the piston, always press it down to the lower stop (parking position). If the piston is not pressed down to the lower stop, unintentional media leaks may occur.

## 5.2 Accessories

### 5.2.1 Flexible discharge tube with recirculation valve

For serial dispensing, the flexible discharge tubing for the bottle-top dispensers Dispensette® S and Dispensette® S Organic can also be used (Accessories/spare parts, p. 82).

The specified accuracy and coefficient of variation of the instrument are only obtained for volumes > 2 ml and by gently approaching the upper and lower stops. The coil of the tubing can be stretched to a maximum length of 800 mm. Before use, the tube must lie in regular loops and must not be twisted. The operating exclusions for the respective instrument in use apply.

### Assembly

#### ⚠ WARNING

#### Use only undamaged tubing


The tubing may not show any signs of damage (e.g. kinks). This must be checked carefully before each use.

- If aggressive liquids are to be dispensed, we recommend using a face shield in addition to the usual safety precautions.
- The bottle must be secured with a bottle holder.
- In order to avoid splashing reagent, always hold the dosing tube firmly, and insert it into the holder provided after use.
- Rinse the tubing to clean it.
- Do not disassemble.

#### Prerequisite:

- If the instrument was in use, it must be cleaned before installing the flexible discharge tube (Cleaning, p. 70).
- a. Instruments with a recirculating valve should be set to 'Recirculate', and the valve lever pulled upwards to remove it.
- b. Slide the discharge tube housing all the way up, then pull it forward with gentle up and down motions.

- c. Push the flexible discharge tube holder from the bottom of the valve block and tighten it. In doing so, the instrument must not be mounted to a bottle. Install the receiver tube.
- d. Press the plug of the recirculation valve downwards.
- e. Slide the flexible discharge tube housing into the valve block up to the stop.
- f. Push the housing all the way down.
- g. Put on the valve toggle made for the discharge valve and press it in firmly. In doing so, pay attention to the color coding and label (see assembly instructions 'Flexible discharge tubing for Dispensette® S').

#### NOTICE

Use a bottle stand (Accessories/spare parts, p. 82).

## 5.2.2 Drying tube

Use of a drying tube, filled with a suitable absorbent (purchased separately), might be necessary for moisture- and CO<sub>2</sub>- sensitive media.

(Accessories/spare parts, p. 82)

### Assembly

- a. Use the assembly tool to unscrew the air vent cap.
- b. Screw in the filled drying tube.
- c. Place the PTFE sealing ring on the bottle thread or the screwed-on adapter and screw the instrument onto the bottle.

#### NOTICE

If necessary, seal the threads of the drying tube, the bottle and/or the bottle adapter with PTFE tape.

## 5.2.3 Sealing ring for valve block

For highly volatile media we recommend to seal the connection from valve block to bottle with the PTFE sealing ring and PTFE tape (Accessories/spare parts, p. 82).

### Assembly

- a.** Place the PTFE sealing ring on the bottle thread or the screwed-on adapter and screw the instrument onto the bottle.

## 5.2.4 Air vent cap for micro filter with Luer-cone

For sterile media we recommend the air vent cap with Luer-cone to attach a micro filter. This provides increased protection against contamination by displacement air (Accessories/spare parts, p. 82).

### Assembly

- a.** Unscrew the air vent cap (see 'Assembly Drying tube', Fig. a.).
- b.** Screw in the air vent cap with a Luer cone.
- c.** Place the PTFE sealing ring on the bottle thread or the screwed-on adapter and screw the instrument onto the bottle.
- d.** Insert a commercially available sterile filter into the Luer cone.

## 6 Error limits

Error limits in relation to the nominal capacity (= maximum volume) indicated on the instrument at equal temperature (20 °C/68 °F) of instrument, ambient environment and distilled water. Testing takes place according DIN EN ISO 8655-6 with a completely filled instrument and with uniform and smooth dispensing.

### Error limits

Nominal volume ml	A* ± %	µl	CV* ≤%	µl
1	0.5	5	0.1	1
2	0.5	10	0.1	2
5	0.5	25	0.1	5
10	0.5	50	0.1	10
25	0.5	125	0.1	25
50	0.5	250	0.1	50
100	0.5	500	0.1	100

$$A_T = \frac{V_N}{V_T} \cdot A_N$$

#### Partial volume

The percentage values for A and CV are relative to the nominal volume ( $V_N$ ) and must be converted for partial volumes ( $V_T$ ).

i.e.	Volume	A* ± %	µl	CV* ≤%	µl
$V_N$	25.0	0.5	125	0.1	25
$V_p = 50\% N$	12.5	1.0	125	0.2	25
$V_p = 10\% N$	2.5	5.0	125	1.0	25

\*A = Accuracy, CV = Coefficient of Variation

### NOTICE

The error limits in DIN EN ISO 8655-5 are significantly lower. The maximum error for a single measurement can be approximated from the sum of error limits  $EL = A + 2 \times CV$  (for 25 ml:  $125 \mu\text{l} + 2 \times \mu\text{l} = 175 \mu\text{l}$ ).

## 7 Checking the Volume (Calibration)

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3–12 months. This time frame should be adjusted to correspond with individual requirements.

Gravimetric volume testing according to DIN EN ISO 8655-6 (for measurement conditions, see Error limits, p. 66) is performed as follows:

### 1. Prepare the instrument

Clean the instrument (Cleaning, p. 70), fill with distilled water, and prime carefully.

### 2. Check the volume

- a. 10 dispensing operations with distilled H<sub>2</sub>O in 3 volume ranges (100%, 50%, 10%) are recommended
- b. For discharge, depress the piston slowly and steadily without force until the lower stop
- c. Wipe off the tip of the dispensing tube.
- d. Weigh the dispensed quantity on an analytical balance. (Please refer to the user manual of the scale manufacturer.)
- e. Calculate the dispensed volume. The Z factor takes account of the temperature and air buoyancy.

### Calculation (for nominal volume)

$x_i$  = weighing results

$n$  = number of weighings

$V_0$  = nominal volume

$Z$  = Correction factor (e.g. 1.0029 µl/mg at 20°C, 1013 hPa)

Mean:

$$\bar{x} = \frac{\sum x_i}{n}$$

Mean volume:

$$\bar{V} = \bar{x} * Z$$

Accuracy\*:

$$A\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Coefficient of variation\*:

$$CV\% = \frac{100 s}{\bar{V}}$$

Standard deviation\*:

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

\*) Accuracy and coefficient of variation are calculated according to the formulas of statistical quality control.

#### NOTICE

Test Instructions (SOPs) are available for download from [www.brand.de](http://www.brand.de).

## 8 Adjustment

After a long period of usage, an adjustment of the instrument might be necessary.

- a. Calibrate for example at nominal volume (Checking the Volume (Calibration)).
- b. Calculate mean volume (actual value) (Checking the Volume (Calibration)).
- c. Adjust the instrument (to the calculated mean volume).
- d. After the adjustment, further calibration is necessary to confirm the appropriate adjustment.


### Example

The gravimetric check gives an actual value of 9.90 ml for a 10 ml instrument set for a nominal volume of 10.00 ml.


### 8.1 Digital model

- a. Slide the latch to the left and remove the front housing shell.

- b. Pull out the safety lock. The adjustment cover will then come off. Discard the adjustment cover.

- c. Pull the red knob to disengage the gears. Set the display to actual delivered volume (e.g., 9.90 ml).

- d. Press in the red knob and then the safety lock again.

- e. Replace housing and slide the latch to the right.
  - ⇒ Adjustment is completed.
  - ⇒ Alteration of the factory setting is indicated by a red recalibration flag (circle).

## Adjustment range

Nominal volume [ml]	Digital max. +/- [ $\mu$ l]
1	12
2	24
5	60
10	120
25	300
50	600
100	–

## 8.2 Analog-adjustable model

- a.** Insert the pin of the assembly tool into the cover plate (Pos. 1), and break it off with a rotating motion. Discard the adjustment cover.

- b.** Insert the pin of the assembly tool into the adjustment screw and rotate to the left in order to increase the dispensing volume, or rotate to the right to decrease the dispensing volume (e.g. for an actual value of 9.97 ml, rotate approx. 1/2 turn to the left).

- c.** Adjustment is completed.  
 ⇒ The change in the adjustment is indicated by a red disk (circle in figure).

## Adjustment range

Nominal volume	Analog-adjustable/Fixed-volume max. +/- [ $\mu$ l]	One rotation corresponds to [ $\mu$ l]
1 ml	6	~ 8
2 ml	12	~ 16
5 ml	30	~ 40
10 ml	60	~ 80
25 ml	150	~ 130
50 ml	300	~ 265
100 ml	600	~ 400

## 9 Cleaning

### ⚠ WARNING

#### Components filled with reagent

The cylinder, valves, telescopic filling tube and discharge tube contain reagent!

- > Never remove the discharge tube while the dispensing cylinder is filled.
- > Point the valves and tube openings away from your body.
- > Wear protective clothing, eye protection and protective gloves!

The instrument must be cleaned in the following situations to assure correct operation:

- immediately if the piston becomes stiff
- prior to a reagent change
- prior to long term storage
- prior to dismantling the instrument
- before autoclaving
- before valve replacement
- regularly when using liquids which form deposits (e.g., crystallizing liquids)
- regularly if liquids accumulate in the screw cap
- regularly if using hygroscopic media (e.g., concentrated sulphuric acid)

For the proper cleaning and removal of any deposits in the parts through which liquid flows, the piston must always be completely removed from the cylinder after rinsing with a suitable cleaning solution. If necessary, the parts can also be cleaned in an ultrasonic bath.

## 9.1 Cleaning, Analog-adjustable model, Fixed-volume model

### 1. Completely emptying the instrument

- a. Screw the instrument onto an empty bottle and empty it completely by dispensing. If the instrument is equipped with a recirculation valve, it must be emptied in both the 'dispense' and 'recirculate' settings.

### 2. Rinsing the instrument

- a. Screw the instrument onto a bottle filled with a suitable cleaning agent (e.g. deionized water) and rinse the instrument several times by completely filling and emptying it.

### 3. Removing pistons


#### NOTICE

#### Avoid swapping pistons

The pistons and cylinders are individually matched, and should not be swapped with pistons from other instruments!

- a. Hold the housing securely and unscrew the piston seat completely by turning it to the left.

- b. Pull out the piston carefully.

#### 4. Clean the piston and cylinder, and refit the device

##### NOTICE

##### Special feature of Dispensette® S Organic

With the Dispensette® S Organic, always insert the piston vertically into the cylinder while rotating it.

- a. Clean the piston and cylinder. Carefully remove any deposits from the upper edge of the dosing cylinder (e.g., with a bottle brush and thin soapy water).
- b. Rinse the piston and cylinder with deionized water, and dry carefully.
- c. Push the piston fully into the cylinder, and reassemble the device. The piston can be inserted into the cylinder only if the dosing cannula is attached.

## 9.2 Cleaning Digital model

### 1. Completely emptying the instrument

- a. Screw the instrument onto an empty bottle and empty it completely by dispensing. If the instrument is equipped with a recirculation valve, it must be emptied in both the 'dispense' and 'recirculate' settings.

## 2. Rinsing the instrument


- a. Screw the instrument onto a bottle filled with a suitable cleaning agent (e.g. deionized water) and rinse the instrument several times by completely filling and emptying it.

## 3. Removing pistons

### NOTICE

#### Avoid swapping pistons

The pistons and cylinders are individually matched, and should not be swapped with pistons from other instruments!

- a. Carry out assembly and disassembly at the maximum volume setting only. To do this, turn the volume adjustment wheel as far as it will go in + direction (= maximum volume).
- b. Slide the latch to the left.
- c. Remove the front housing shell.
- d. Use the assembly tool to loosen the piston mounting nut and pull the rear housing section backwards.
- e. Then carefully pull the piston out of the cylinder.

## 4. Clean the piston and cylinder, and refit the device

### NOTICE

#### Special feature of Dispensette® S Organic

With the Dispensette® S Organic, always insert the piston vertically into the cylinder while rotating it.



- a. Clean the piston and cylinder. Carefully remove any deposits from the upper edge of the dosing cylinder (e.g., with a bottle brush and thin soapy water).
- b. Rinse the piston and cylinder with deionized water, and dry carefully.
- c. Push the piston fully into the cylinder, and reassemble the device. The piston can be inserted into the cylinder only if the dosing cannula is attached.

## NOTICE

### Cleaning Digital model

The stop segment must engage underneath the cylinder's stop ring. When tightening the piston mounting nut using the assembly tool, use your thumbs to press the piston/cylinder unit back towards the back side of the housing.

## 9.3 Replacing discharge tube/valves

### NOTICE

A functional test must always be carried out after exchanging components.

### 9.3.1 Replacing discharge tube

- a. Instruments with a recirculating valve should be set to 'Recirculate', and the valve lever pulled upwards to remove it.

- b. Slide the discharge tube housing all the way up, then pull it forward with gentle up and down motions.
- c. Hold coupling piece of the new discharge tube and pull housing up. Push housing into the valve block until it meets the stop.

- d. Slide the discharge tube housing all the way down.
- e. For instruments with a recirculation valve, pull up the valve lever to the 'Recirculate' position, and press it in tightly.

## 9.3.2 Replacing valves

### NOTICE

Always install the valve intended for the particular instrument model and size! (See ordering data, Ordering Information, p. 79) Identical filling valves but different discharge valves are used for Dispensette® S and Dispensette® S Organic. For easy identification, the discharge valves for the Dispensette® S Organic are marked with 'ORG'.

### 9.3.2.1 Replacing the discharge valve


#### NOTICE

The discharge valve has a safety ball that closes automatically when the discharge tube is removed. Attaching the discharge tube reopens the safety ball.

- a. After disassembling the discharge tube (see Replacing discharge tube, p. 73), use the assembly tool to unscrew the discharge valve.
- b. Firmly screw in the new discharge valve first by hand and then tighten with the assembly tool. The thread should not be visible.

### 9.3.2.2 Replacing the filling valve

- a. Pull out the recirculation tube and the telescopic filling tube.

- b. Use the assembly tool to unscrew the filling valve.
- c. Screw in the new filling valve first by hand and then tighten it with the assembly tool.

### 9.3.2.3 Dislodging a stick valve ball


If the instrument does not fill up, and if some elastic resistance is evident when the piston is pulled upward, then it is possible that the valve ball is stuck.

In this case, dislodge the valve ball using light pressure; for example, with a 200 µl plastic pipette tip (see the figure at the side).

# 10 Autoclaving

The device can be autoclaved at 121°C (250°F), 2 bar, and a holding time of at least 15 minutes according to DIN EN 285.

- a. The instrument must be carefully cleaned prior to autoclaving (Cleaning, p. 70).
- b. Open the closure cap on the discharge tube, and for instruments with a recirculation valve, set the valve to 'Dispense'.
- c. Check that the filling valve is securely seated (Fig. 1a). In the Digital model, additionally check that the piston mounting is securely seated (Fig. 1b).
- d. To ensure unhindered access for the steam and to prevent the valve ball in the filling valve from possibly becoming stuck, hold the instrument with the discharge piston pressed vertically downward, and gently tap against the casing with your hand (Fig. 2). Then lay it horizontally in the autoclave. Be sure to avoid the instrument coming into contact with metal surfaces in the autoclave!

## NOTICE

Do not reassemble the instrument until it has cooled down to room temperature (Cooling time approx. 2 hours). After every autoclaving, inspect all parts for deformities or damage. If necessary, replace them. The effectiveness of autoclaving must be verified by the user.

# 11 Troubleshooting

Problem	Possible cause	Corrective action
Piston sluggish or stuck	Crystal deposits Contaminants	Stop dispensing immediately. Loosen piston with circular motion, but do not disassemble. Perform cleaning (Cleaning, p. 70).
Filling not possible	Volume adjusted to minimum setting	Set the desired volume (Dispensing, p. 62).
	Filling valve stuck	Unscrew the filling valve from the valve block and clean it. If the valve ball is stuck, use a 200 µl plastic tip to dislodge it (Dislodging a stick valve ball, p. 74). If necessary, replace the filling valve.
Dispensing not possible	Discharge valve stuck	Unscrew the discharge valve from the valve block, clean it and replace if necessary. If the valve ball is stuck, use a 200 µl plastic tip to dislodge it.
Discharge tube (with or without recirculation valve) cannot be mounted sufficiently	Discharge valve is not screwed in deeply enough	Tighten the discharge valve with the assembly tool until it meets the stop so that the threads are no longer visible.
Air bubbles in the instrument	Reagent with high vapor pressure has been drawn in too quickly	Slowly draw in reagent.
	Valve screw connections loose	Tighten the valves firmly with the assembly tool.
	The instrument has not been primed	Prime the instrument (Priming, p. 60).
	Filling tube is loose or damaged	Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube.
	Valves not firmly connected or damaged	Perform cleaning (Cleaning, p. 70). Tighten the valves using the assembly tool.
Dispensed volume too low	Filling tube is loose or damaged	Perform cleaning (Cleaning, p. 70). Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube (Replacing discharge tube, p. 73).
	Filling valve is loose or damaged	Perform cleaning (Cleaning, p. 70). Tighten the valves using the assembly tool. If necessary, replace filling valves.

<b>Problem</b>	<b>Possible cause</b>	<b>Corrective action</b>
Leaking liquid between instrument and bottle	Recirculation tube not connected	Connect recirculation tube (First Steps, p. 59, Fig. 3).
	Volatile reagent dispensed without sealing ring	Install sealing ring for valve block (First Steps, p. 59).
Liquid discharge at the upper edge of the dispensing cylinder	Piston not in lowest position after dispensing	After use, always depress the piston to the lower stop (put it in the "parking position"). (see Dispensing, p. 62)
	Hygroscopic dispensing medium leads to a leakage of medium	Regularly wipe off any accumulation of liquid on the upper edge of the dispensing cylinder

## 12 Product markings

Symbol or number	Meaning
	General warning sign
	Note user manual
	Note eye protection
	Use hand protection
	Use protective clothing
XXXXXXXX	Serial number
	The device is marked in accordance with the German Weights and Measures Act and the Weights and Measures Ordinance. Character sequence DE-M (DE for Germany), framed by a rectangle, as well as the two last digits of the year the marking was added.
www.brand.de/ip	Patent information

# 13 Ordering Information

## Dispensette® S, Digital

Volume ml	Graduation ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
0.1 – 1	0.005	<a href="#">4600310</a>	<a href="#">4600311</a>
0.2 – 2	0.01	<a href="#">4600320</a>	<a href="#">4600321</a>
0.5 – 5	0.02	<a href="#">4600330</a>	<a href="#">4600331</a>
1 – 10	0.05	<a href="#">4600340</a>	<a href="#">4600341</a>
2.5 - 25	0.1	<a href="#">4600350</a>	<a href="#">4600351</a>
5 – 50	0.2	<a href="#">4600360</a>	<a href="#">4600361</a>

## Dispensette® S, Analog-adjustable

Volume ml	Graduation ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
0.1 - 1	0.02	<a href="#">4600100</a>	<a href="#">4600101</a>
0.2 - 2	0.05	<a href="#">4600120</a>	<a href="#">4600121</a>
0.5 - 5	0.1	<a href="#">4600130</a>	<a href="#">4600131</a>
1 - 10	0.2	<a href="#">4600140</a>	<a href="#">4600141</a>
2.5 - 25	0.5	<a href="#">4600150</a>	<a href="#">4600151</a>
5 - 50	1.0	<a href="#">4600160</a>	<a href="#">4600161</a>
10 - 100	1.0	<a href="#">4600170</a>	<a href="#">4600171</a>

## Dispensette® S, Fixed-volume

Volume ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
1	<a href="#">4600210</a>	<a href="#">4600211</a>
2	<a href="#">4600220</a>	<a href="#">4600221</a>
5	<a href="#">4600230</a>	<a href="#">4600231</a>
10	<a href="#">4600240</a>	<a href="#">4600241</a>
Fixed volumes by request: 0.5–100 ml (please specify when ordering)	<a href="#">4600290</a>	<a href="#">4600291</a>

**NOTICE**

Scope of delivery, see Scope of delivery, p. 49

**Dispensette® S Organic, Digital**

Volume ml	Graduation ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
0.5 - 5	0.02	<a href="#">4630330</a>	<a href="#">4630331</a>
1 - 10	0.05	<a href="#">4630340</a>	<a href="#">4630341</a>
2.5 - 25	0.1	<a href="#">4630350</a>	<a href="#">4630351</a>
5 - 50	0.2	<a href="#">4630360</a>	<a href="#">4630361</a>

**Dispensette® S Organic, Analog-adjustable**

Volume ml	Graduation ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
0.5 - 5	0.1	<a href="#">4630130</a>	<a href="#">4630131</a>
1 - 10	0.2	<a href="#">4630140</a>	<a href="#">4630141</a>
2.5 - 25	0.5	<a href="#">4630150</a>	<a href="#">4630151</a>
5 - 50	1.0	<a href="#">4630160</a>	<a href="#">4630161</a>
10 - 100	1.0	<a href="#">4630170</a>	<a href="#">4630171</a>

**Dispensette® S Organic, Fixed-volume**

Volume ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
5	<a href="#">4630230</a>	<a href="#">4630231</a>
10	<a href="#">4630240</a>	<a href="#">4630241</a>
Fixed volumes by request: 2-100 ml (please specify when ordering)	<a href="#">4630290</a>	<a href="#">4630291</a>

**NOTICE****Dispensing hydrofluoric acid (HF)**

For dispensing hydrofluoric acid (HF), we recommend the use of the Dispensette® S Trace Analysis bottle-top dispenser with platinum-iridium valve spring (see separate operating manual).

---

# 14 Accessories/spare parts

## Bottle adapter

External thread	for bottle thread/ground joint	Material	Order No.
GL 32	GL 25	PP	<a href="#">704325</a>
GL 32	GL 28/ S 28	PP	<a href="#">704328</a>
GL 32	GL 30	PP	<a href="#">704330</a>
GL 32	GL 45	PP	<a href="#">704345</a>
GL 45	GL 32	PP	<a href="#">704396</a>
GL 45	GL 35	PP	<a href="#">704431</a>
GL 45	GL 38	PP	<a href="#">704397</a>
GL 45	S* 40	PP	<a href="#">704343</a>
GL 45	S* 50	PP	<a href="#">704350</a>
GL 45	S* 54	PP	<a href="#">704430</a>
GL 45	S* 60	PP	<a href="#">704348</a>
GL 32	GL 25	ETFE	<a href="#">704375</a>
GL 32	GL 28/ S 28	ETFE	<a href="#">704378</a>
GL 32	GL 45	ETFE	<a href="#">704395</a>
GL 45	GL 32	ETFE	<a href="#">704398</a>
GL 45	GL 38	ETFE	<a href="#">704399</a>
GL 45	S* 40	PTFE	<a href="#">704391</a>
GL 32	NS 19/26	PP	<a href="#">704419</a>
GL 32	NS 24/29	PP	<a href="#">704424</a>
GL 32	NS 29/32	PP	<a href="#">704429</a>

\* buttress thread

### NOTICE

#### Template for adapter

An adapter template is available for download at [www.brand.de](http://www.brand.de) in the Service section. This can be used to determine thread and adapter sizes. This adapter template for the respective product can also be found in the on-line shop [store.brand.de](http://store.brand.de) in the “Selection Charts” section.

## Discharge tubes without recirculation valve

### Discharge tube without recirculation valve for Dispensette® S

Pkg unit Pack of 1.

Nominal volume ml	Model	Length (mm)	Order no.
1, 2, 5, 10	Fine tip	108	<a href="#">708002</a>
5, 10	Standard	108	<a href="#">708005</a>
25, 50, 100	Fine tip	135	<a href="#">708006</a>
25, 50, 100	Standard	135	<a href="#">708008</a>

## Discharge tubes without recirculation valve for Dispensette® S Organic

Pkg unit Pack of 1.

Nominal volume ml	Model	Length (mm)	Order no.
5, 10	Fine tip	108	<a href="#">708012</a>
5, 10	Standard	108	<a href="#">708014</a>
25, 50, 100	Fine tip	135	<a href="#">708016</a>
25, 50, 100	Standard	135	<a href="#">708019</a>

## Discharge tubes with recirculation valve

### Discharge tubes with recirculation valve for Dispensette® S

Pkg unit Pack of 1.

Nominal volume ml	Model	Length (mm)	Order no.
1, 2, 5, 10	Fine tip	108	<a href="#">708102</a>
5, 10	Standard	108	<a href="#">708104</a>
25, 50, 100	Fine tip	135	<a href="#">708106</a>
25, 50, 100	Standard	135	<a href="#">708109</a>

## Discharge tubes with recirculation valve for Dispensette® S Organic

Pkg unit Pack of 1.

Nominal volume ml	Model	Length (mm)	Order no.
5, 10	Fine tip	108	<a href="#">708112</a>
5, 10	Standard	108	<a href="#">708114</a>
25, 50, 100	Fine tip	135	<a href="#">708116</a>
25, 50, 100	Standard	135	<a href="#">708119</a>

## Flexible discharge tubing with recirculation valve for Dispensette® S and Dispensette® S Organic

PTFE, spiraled, approx. 800 mm long with safety handle.

Pkg unit Pack of 1.

Not suitable for hydrofluoric acid (HF)

Nominal volume ml	Discharge tubing outer diameter	Discharge tubing inner diameter	Order no.
1, 2, 5, 10	3	2	<a href="#">708132</a>
25, 50, 100	4.5	3	<a href="#">708134</a>

## Discharge valve Dispensette® S

PFA/Boro 3.3/ceramic/platinum-iridium.

No valve marking.

Pkg unit Pack of 1.

for nominal volume ml	Order no.
1, 2*	<a href="#">6749</a>
5, 10	<a href="#">6727</a>
25, 50, 100	<a href="#">6728</a>

\* with valve marking '1 + 2'

## Dispensette® S Organic discharge valve

PFA/Boro 3.3/ceramic/tantalum.

Valve marked with 'ORG'.

Pkg unit Pack of 1.

for nominal volume ml	Order no.
5, 10	<a href="#">6729</a>
25, 50, 100	<a href="#">6730</a>

## Dispensette® S and Dispensette® S Organic filling valve


PFA/ETFE/Boro 3.3/ceramic.

No valve marking.

Pkg unit Pack of 1.

for nominal volume ml	Order no.
1, 2, 5, 10	<a href="#">6734</a>
25, 50, 100	<a href="#">6735</a>

## Telescopic filling tubes for Dispensette® S and Dispensette® S Organic


FEP.

Individually adjustable lengths.

Packaging unit Pack of 1.

for nominal volume ml	Outer diameter mm	Length (mm)	Order no.
1, 2, 5, 10	6	70-140	<a href="#">708210</a>
		125-240	<a href="#">708212</a>
		195-350	<a href="#">708214</a>
		250-480	<a href="#">708216</a>
25, 50, 100	7.6	170-330	<a href="#">708218</a>
		250-480	<a href="#">708220</a>

## Venting stopper for microfilter with Luer cone

Description	Pkg unit	Order no.
Venting stopper for microfilter with Luer cone. PP. Air vent cap and PTFE-sealing ring.	1 pc.	<a href="#">704495</a>

## Recirculation tube

Description	Pkg unit	Order no.
Recirculation tube. FEP	1 pc.	<a href="#">6747</a>

## Adjustment/installation tool

Description	Pkg unit	Order no.
Adjustment/installation tool	1 pc.	<a href="#">6748</a>

## Bottle stand


Description	Pkg unit	Order no.
Bottle stand PP. Support rod, 325 mm, base plate 220 x 160 mm.	Pack of 1	<a href="#">704275</a>

## Sealing ring for valve block

Description	Pkg unit	Order no.
Sealing ring for valve block. PTFE, for highly volatile media.	1 pc.	<a href="#">704486</a>

## Drying tube

Description	Pkg unit	Order no.
Drying tube, including PTFE sealing ring (without drying agent)	Pack of 1	<a href="#">707930</a>


## Screw cap with fastener

Packaging unit Pack of 1.

Description	Nominal volume ml	Order No.
PP, red, for Dispensette® S	1, 2, 5, 10	<a href="#">706018</a>
	25, 50, 100	<a href="#">706019</a>
PP, yellow, for Dispensette® S Organic	5, 10	<a href="#">706025</a>
	25, 50, 100	<a href="#">706027</a>
ETFE, for Dispensette® S and Dispensette® S Organic	1, 2, 5, 10	<a href="#">706029</a>
PTFE, for Dispensette® S and Dispensette® S Organic	25, 50, 100	<a href="#">706031</a>

## Dispensette® drum extraction system

Pkg unit Pack of 1.

Description	Order no.
Dispensette® drum extraction system	<a href="#">704261</a>
Tripod connection for wall mount, stainless steel, Dispensette® extraction system	<a href="#">704268</a>
Table/rack clamp for wall mount for Dispensette® extraction system	<a href="#">704272</a>

# 15 Repairs

## 15.1 Sending for repair

### NOTICE

Transporting of hazardous materials without a permit is a violation of federal law.

### Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

### Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from [www.brand.de](http://www.brand.de).

### Outside USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

### Contact addresses

#### Germany:

BRAND GMBH + CO KG  
Otto-Schott-Straße 25  
97877 Wertheim (Germany)  
T +49 9342 808 0  
F +49 9342 808 98000  
[info@brand.de](mailto:info@brand.de)  
[www.brand.de](http://www.brand.de)

#### USA and Canada:

BrandTech® Scientific, Inc.  
11 Bokum Road  
Essex, CT 06426-1506 (USA)  
T +1-860-767 2562  
F +1 - 860 - 767 2563  
[info@brandtech.com](mailto:info@brandtech.com)  
[www.brandtech.com](http://www.brandtech.com)

#### India:

BRAND Scientific Equipment Pvt. Ltd.  
303, 3rd Floor, 'C' Wing, Delphi  
Hiranandani Business Park,  
Powai  
Mumbai-400 076 (India)  
T +91 22 42957790  
F +91 22 42957791

#### China:

BRAND (Shanghai) Trading Co., Ltd.  
Rm 201-202, North Tower,  
No. 199 Kaibin Rd, Xuhui District, Shanghai  
Shanghai 200030 (P.R. China)  
T +86 21 6422 2318  
F +86 21 6422 2268  
[info@brand.com.cn](mailto:info@brand.com.cn)

info@brand.co.in  
www.brand.co.in

www.brand.cn.com

## 16 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring devices. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the device. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive media.

The complete SOP for testing can be downloaded from [www.brand.de](http://www.brand.de) or [www.brandtech.com](http://www.brandtech.com).

BRAND also offers you the option of having your devices calibrated through our factory calibration service or through our accredited calibration laboratory. Just send us the devices to be calibrated, indicating the type of calibration you would like. You will get your devices back in a few days. A detailed calibration report (factory calibration) or an accredited calibration certificate in accordance with DIN EN ISO/IEC 17025 is enclosed with each device. More information can be obtained from your retailer or directly from BRAND. The order document is available for download at [www.brand.de](http://www.brand.de) (Service & Support).

### For customers outside Germany

If you would like to use our calibration service, please contact one of our service partners in your region. Our service partners can forward your devices to BRAND for factory calibration, if required.

## 17 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

### **USA and Canada:**

Find more warranty information on [www.brandtech.com](http://www.brandtech.com).

## 18 Disposal

Before disposal, observe the relevant national disposal regulations, and ensure that the product is disposed of properly.

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