Liquid Handling Station

Pipettierroboter | Pipetting robot

Impressum

BRAND GMBH + CO KG
Otto-Schott-Str. 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

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i.A. Siegfried Ott
Regulatory Affairs

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UK Declaration of Conformity

Trade name of the device: Liquid Handling Station incl. Liquid Ends
 Device variations: Liquid Handling Station with FlowBox
 Cat.-No.: 709400; 709402
 Manufacturer: BRAND GMBH + CO KG
 Otto-Schott-Str. 25
 97877 Wertheim · Germany

We, as the manufacturer of the above described device(s) take sole responsibility for and hereby declare that the described device(s) meet(s) the provisions of the following Regulation(s)/Directive(s):	Applied harmonized standards:
The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012	EN IEC 63000:2018
Electromagnetic Compatibility Regulations 2016 (SI 2016 No. 1091, as amended by SI 2019 No. 696)	EN 61326-1:2013
Electrical Equipment (Safety) Regulations 2016 (SI 2016 No. 1101, as amended by SI 2019 No. 696)	EN 61010-1:2010, EN 61010-1:2010/ A1:2019/AC:2019-04, EN 61010-1:2010/ A1:2019
Other applied standards:	Scope:
EN IEC 61010-2-081:2020	Electrical Equipment (Safety) Regulations 2016

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S. Ott

i.A. Siegfried Ott
 Regulatory Affairs

This document declares the accordance with the named harmonized regulations, but does NOT assure specific properties.

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Table of Contents

1	Devices and components	55
1.1	Scope of supply.....	55
1.2	Safety	55
1.2.1	Intended use.....	55
1.2.2	Safety instructions	58
1.2.3	Transport and storage	61
1.3	Liquid Handling Station and accessories	62
1.3.1	Liquid Handling Station.....	62
1.3.2	Liquid Handling Station flow.....	63
1.3.3	Liquid End.....	64
1.3.4	Teleshake 95.....	64
1.4	Technical data	66
1.5	Accuracy Table.....	68
1.6	Commissioning	68
1.6.1	Installation	68
1.6.2	Connection	70
1.6.3	Connecting the Liquid Handling Station	70
1.6.4	Connecting the Liquid Handling Station flow	71
1.6.5	Connecting the shaking and temperature control module Teleshake 95	71
1.7	Operation	72
1.7.1	Liquid Handling Station.....	72
1.7.2	Liquid Handling Station flow.....	73
1.7.3	Liquid End.....	74
1.7.4	Shaking and temperature control module Teleshake 95	74
1.8	Maintenance/Cleaning	75
1.8.1	Device	75
1.8.2	Liquid Ends.....	76
1.9	Ordering Information	84
1.9.1	Liquid Handling Station and accessories	84
1.9.2	Spare parts	90
2	Glossary.....	93

1 Devices and components

1.1 Scope of supply

Scope of supply for Liquid Handling Station

The scope of supply includes:

- BRAND Liquid Handling Station
- Motor control unit
- Operating software
- Operating manual
- Power cable
- USB cable
- Documentation and on-site training

Scope of supply for Liquid Handling Station flow

The scope of supply includes:

- BRAND Liquid Handling Station flow
- Protective screen for inserting in the Liquid Handling Station flow
- Motor control unit
- Operating software
- Operating manual
- 2 x power cable
- 2 x USB cable
- 24 V DC universal power adapter for FlowBox
- Documentation and on-site training

1.2 Safety

1.2.1 Intended use

1.2.1.1 Liquid Handling Station

The BRAND Liquid Handling Station is a pipetting robot for use in routine laboratories with small and medium sample throughput in research, development, and production, especially in fields of application like PCR, qPCR, ELISA, enzyme tests, etc. The device is typically used to pipette aqueous media such as buffer solutions (phosphate buffer, Tris/HCl buffer, etc.), protein solutions (BSA solutions, enzyme solutions, PCR master mixes), and aqueous media and samples.

The device is used for the automatic, precise transfer of liquids. For the liquid transfer, autoclavable 1-channel and 8-channel pipetting modules (Liquid Ends) are available that must be changed manually. 7 SLAS work stations (P2-P8) can be freely equipped on the work table. An additional predefined work station (P1) must be assigned with the waste container.

1.2.1.2 Liquid Handling Station flow

⚠ DANGER

Use of harmful substances, infectious liquids or flammable media

When working with the FlowBox, a continuous stream of air is blown out of the pipetting robot into the environment. This air may contain harmful or infectious aerosols.

- The FlowBox may not be switched on if you work with biohazardous, radioactive, toxic or otherwise hazardous substances!
- Do not use any media in the work area and in the area around the FlowBox that may generate strong heat in a reaction or form a potentially explosive atmosphere. Examples: Acid-base reactions, carbon disulfide.
- The device must be disinfected prior to using the FlowBox according to the instructions in the chapter Maintenance/Cleaning. If you have any questions, please contact the manufacturer and/or your safety officer.

NOTICE

Turbulence

In order to reduce turbulence, open the front door slowly and as rarely as possible. Also, the protective screen must be inserted.

NOTICE

Regularly exchanging the HEPA filter

The HEPA filter of the FlowBox gets dirty with time. It is therefore recommended to have the filter changed by a service employee at least once per year.

The BRAND Liquid Handling Station flow is a pipetting robot with a permanently installed filter unit, the so-called FlowBox, which purifies the air in the work area of the Liquid Handling Station flow. For this purpose, ambient air is aspirated, filtered, and blown outwards in a laminar flow from back to front across the work table through small openings in the front door. Particles that entered the work area of the Liquid Handling Station flow before switching on the FlowBox are flushed out of the work area through the air stream and the openings in the front door. The filter of the FlowBox is therefore used to reduce the number of particles in the work area. If the front door of the operationally ready pipetting robot is opened, the fan speed and thus the air stream directed outwards from the FlowBox are increased in order to reduce the risk of sample contamination.

1.2.1.3 Usage limits and restrictions

Limits	Limit
Operating temperature	+ 15 ... + 35 °C (device and reagents)
Vapor pressure	500 mbar
Viscosity	260 mPa s

The viscosity of liquids as well as wetting liquids compromise the accuracy, as do liquids whose temperature deviates from room temperature by more than +/- 5 °C.

NOTICE

Risk of corrosion

The device can be damaged by media.

- Do not use the device to pipette liquids that attack polypropylene, PMMA (side and front panel), POM or aluminum (labware adapters).
- Avoid aggressive vapors.
- Avoid strongly acidic or alkaline solutions.

Protect the device from the penetration of liquids. Do not place any vessels containing liquids on or next to the device.

This device can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

In case of doubt, contact the manufacturer.

The safety instructions listed in the following chapter must be strictly observed. Claims of any type are excluded for damage resulting from improper use.

1.2.2 Safety instructions

1.2.2.1 Danger from electric current

DANGER

Danger to life from electrical voltage

Life-threatening electrical voltages are generated within the device. Touching live parts poses a danger to life due to electric shock.

- Commissioning may only be carried out if the device has been properly installed / repaired.
- In case of danger, immediately disconnect the device from the mains voltage; to do this, pull the plug or disconnect the cable.
- Ensure at all times that the housing and the covers are intact and closed.
- Only connect the device to sockets with a protective contact.
- Before connecting the device, check whether the mains voltage agrees with the admissible operating voltage. An incorrect mains voltage may result in the destruction of the device.
- Before switching on the device, ensure that the power cable and connections are undamaged.
- The covers of the device may only be removed by qualified persons who are authorized to do so by the manufacturer! Prior to disassembly, switch off the device and pull out the power plug.
- Keep humidity away from live parts in order to prevent short circuits.

1.2.2.2 Personnel requirements

WARNING

Personnel qualification

Unqualified personnel are prohibited from working with the device.

- The device may only be operated by qualified persons, i.e. only by qualified users trained by the manufacturer or specialist dealer.
- In order to prevent injury and material damage, keep unauthorized persons away from the work area of the device.
- Each user must have read the operating manual carefully prior to using the device.

⚠ CAUTION**Repairs and Service**

All repairs and service work may only be carried out by skilled and qualified personnel.

- Repairs and service work may only be performed only by qualified personnel who are authorized by the manufacturer.
- Use only original spare parts.
- Do not make any technical modifications to the device.

1.2.2.3 Danger due to failure to comply with the intended purpose and improper use

⚠ WARNING**Observance of usage limits**

Do not use the device beyond its intended use.

- The device may only be used for its intended purpose. Explosive and flammable substances may not be used.
- The manufacturer assumes no liability for personal injury or material damage resulting from the non-observance of the intended purpose.
- Contact the manufacturer in case of doubt.

⚠ WARNING**Improper handling**

- ✓ Handle the device according to the safety instructions and operating instructions.
- Never use force.
- Improper handling (short circuit, mechanical destruction, overheating, etc.) may potentially result in a fire or the explosion of the device.
- No liability is assumed for damage resulting for improper handling.

1.2.2.4 Locally applicable safety regulations

⚠ WARNING**Compliance with locally applicable regulations**

In addition to the safety regulations for using the device, the locally applicable regulations must be complied with.

- This applies especially for occupational safety and accident prevention (e.g. protective clothing, safety gloves, and eye protection) and hygiene regulations.
- Furthermore, the regulations of your laboratory for the handling of the substances used are also authoritative.

1.2.2.5 Danger from harmful substances or infectious media

DANGER

Use of harmful substances or infectious liquids

When handling infectious aqueous solutions, the national regulations and the biosafety level of the laboratory must be considered. The Liquid Handling Station can be used for handling germs and biological material of **risk class I**. Germs and biological material of **higher risk classes** may not be used.

- The operator bears responsibility for the selection of the media used and their safe handling.
- This especially applies to radioactive, infectious, toxic, aggressive, flammable, and otherwise hazardous substances.
- When working with infectious or hazardous samples, the standard laboratory rules and precautions must be adhered to.
- Observe the material safety data sheets and application notes of the manufacturer.
- Dispose of the content of the waste container according to the applicable regulations. Particular caution must be exercised with flammable and infectious liquids.
- If the device was contaminated by hazardous substances, it must be cleaned and decontaminated.
- No potentially toxic or harmful gases or substances are released by the device itself.

1.2.2.6 Explosion hazard

DANGER

Explosion hazard

In order to prevent the risk of an explosion, comply with the following measures:

- Do not operate the device in potentially explosive atmospheres.
- Do not use explosive or strongly reactive substances with the device.
- Do not store explosive substances in the vicinity of the device.
- Do not use flammable liquids and especially carbon disulfide!

1.2.2.7 Danger from crushing

CAUTION

Crushing hazard - Keep fingers and hands safe!

It is only permitted to reach into the work area once the device has gone into a standstill!

- If the front door is opened during a program, the operation is interrupted once the current work step has been completed.
- Wait until the device is stationary before you reach into the work area!

1.2.2.8 Danger from hot surfaces

⚠ CAUTION

Risk of burns - Keep fingers and hands safe!

- Allow the Teleshake 95 heat and shake unit to cool before you work with it.
- The Teleshake 95 may still be hot after switching off the device.

1.2.2.9 Visual check before beginning work

⚠ CAUTION

Insufficient inspection before beginning work

Before beginning work, make sure that the device is in proper condition.

- If there should be any recognizable damage that threatens safe operation, do not work with the device and arrange for a repair immediately.
- If the device is contaminated by liquids or dirt, remove these contaminants.

1.2.2.10 Malfunctions

⚠ CAUTION

Handling malfunctions

If malfunctions occur during the work process, act immediately.

- Follow the instructions on the screen.
- If the program does not provide any error messages or instructions, stop the process by pressing the on/off switch. Contact the manufacturer in case of doubt.
- An "emergency stop" can be effected by pulling the power plug. Please note that data may be lost in this case.

1.2.3 Transport and storage

⚠ WARNING

Risk of injury and material damage during transport

- At least two people must lift or carry the device.
- When carrying the device, only take hold of the sides of the device from below.
- Only transport the device upright without tipping it.

Prior to transport, remove all loose parts (adapters, vessels, labware, etc.). Do not leave any dangerous substances in the device. If necessary, clean and decontaminate the device prior to transport.

Ambient conditions for transport and storage:

Limits	Limit
Temperature range	- 20 ... + 65 °C
Humidity	max. 95 %

For optimal protection, the device is transported in a special packaging. To preserve this protection during storage, only remove the packaging shortly before installation.


NOTICE

Transport only in original packaging

Always use the original packaging for transport! For this reason, please refrain from damaging the packaging or throwing it away and proceed according to the included instructions when unpacking and packing.

1.3 Liquid Handling Station and accessories

1.3.1 Liquid Handling Station

- 1 Motor control unit
- 2 USB interface
- 3 Power connection


- 5 Support feet, height-adjustable
- 6 On/off switch
- 7 Front door

4 Power button

8 Work Table with positions P1 ... P8


1.3.2 Liquid Handling Station flow

Liquid Handling Station flow

- | | |
|---|--|
| <p>1 On/off switch of the FlowBox</p> <p>2 Motor control unit of the Liquid Handling Station</p> <p>3 Intake filter (side)</p> <p>4 USB interface of the FlowBox</p> <p>5 LED lit: voltage present</p> <p>6 LED lit: Liquid Handling Station flow is ready for operation.</p> <p>7 Power supply IN 24 V DC</p> <p>8 Opening for internal relative pressure measurement. Do not seal off</p> <p>9 USB hub</p> | <p>10 USB interface</p> <p>11 Power connection</p> <p>12 Power button</p> <p>13 On/off switch of the Liquid Handling Station flow</p> <p>14 Protective screen behind front door</p> <p>15 Support feet, height-adjustable</p> <p>16 Front door</p> <p>17 Openings in front door</p> <p>18 Discharge filter in front door</p> |
|---|--|


1.3.3 Liquid End

Ill. 1: Liquid End components

- 1 Decoupling button
- 2 Motor control unit for *Liquid End*
- 3 Liquid End

1.3.4 Teleshake 95

Ill. 2: Teleshake 95

- 1 Adapter plate
- 2 Heat and shake plate
- 3 Sub-D cable
- 4 Centering plate

5 Name plate


1.3.4.1 Single TEC Control

Ill. 3: Single TEC front

- | | |
|--|--|
| <ul style="list-style-type: none"> 1 Switch control unit on/off 3 USB 5 External sensor connection
 7 Slot modules 1 ... 6 for connecting the Tele-shake 95 via Sub-D connector (not drawn in) 9 Touchscreen | <ul style="list-style-type: none"> 2 Power supply 4 ID 6 Color-coding on slot module. The color-coding corresponds with the Sub-D connector that is likewise color-coded. 8 Ventilation openings |
|--|--|

1.3.4.2 Multi TEC Control

Ill. 4: Multi TEC

- | | |
|---|---|
| <p>1 External sensor connection</p> <p>3 USB</p> <p>5 Power supply</p> <p>7 Slot modules 1 ... 6 for connecting the Tele-shake 95 via Sub-D connector</p> <p>9 Touch screen</p> | <p>2 ID</p> <p>4 Switch control unit on/off</p> <p>6 Color-coding on slot module. The color-coding corresponds with the Sub-D connector that is likewise color-coded.</p> <p>8 Ventilation openings</p> |
|---|---|

1.4 Technical data

Liquid Handling Station

Liquid Ends	1-channel Liquid Ends (SC), 8-channel Liquid Ends (MC8)
Volume Range	1-channel Liquid Ends: 1 - 50 µl, 10 - 200 µl, 40 - 1000 µl 8-channel Liquid Ends: 1 - 50 µl, 20 - 300 µl, 40 - 1000 µl
Working positions	8 Working positions: P2 - P8, P1 for waste container
Weight	25 kg
Dimensions	H 530 x W 595 x D 485 mm
Operating temperature	+ 15 °C ... + 35 °C
Transport temperature	- 20 °C ... + 65 °C
Mains voltage	100 ... 240 V, 50 ... 60 Hz
Microfuse	2 x T 2,5A H 250V
Interfaces	1 x USB
Power consumption	max. 150 W
Protection Class	Protection Class I
Housing protection rating	IP20
Safety standards	IEC 61 010-1
EMC compatibility radio interference suppression and interference immunity in accordance with the following standard	DIN EN 61 326-1
Sound level	46 dB

Table 1: Liquid Handling Station

Liquid Handling Station flow

Liquid Ends	1-channel Liquid Ends (SC), 8-channel Liquid Ends (MC8)
-------------	---

Volume Range	1-channel Liquid Ends: 1 - 50 µl, 10 - 200 µl, 40 - 1000 µl 8-channel Liquid Ends: 1 - 50 µl, 20 - 300 µl, 40 - 1000 µl
Working positions	8 Working positions: P2 - P8, P1 for waste container
Weight	approx. 35 kg (pipetting robot plus FlowBox)
Dimensions	H 662 x W 595 x D 515 mm
Operating temperature	+ 15 °C ... + 35 °C
Transport temperature	- 20 °C ... + 65 °C
Mains voltage	100 ... 240 V, 50 ... 60 Hz
Microfuse	2 x T 2,5A H 250V
Interfaces	2 x USB
Power consumption	max. 150 W
Protection Class	Protection Class I
Housing protection rating	IP20
Safety standards	IEC 61 010-1
EMC compatibility	Radio interference suppression and interference immunity in accordance with the following standard: DIN EN 61 326-1
Sound level	46 dB

Flowbox

Sound level	~60 dBA (with front door closed) ~67 dBA (with front door opened)
Volume flow rate	~20 m ³ /h (with front door closed) ~30 m ³ /h (with front door open)
Air change	~260 times per hour (with front door closed) ~330 times per hour (with front door open)
Filter rating	H 14 following EN 1822

Table 2: Flowbox

Flowbox universal power adapter

Input	100-240 V, 50-60 Hz
Output	24V DC, 3,75A DC, 90W max.
Protection Class	Protection Class 1
Safety standard	EN 60950-1

EMC compatibility	EN 55022 class B, EN 61000-4-2,3,4,5,6,8,11
-------------------	---

Table 3: Flowbox universal power adapter

1.5 Accuracy Table

1-Channel Liquid End

Volume range [μl]	Partial volume [μl]	A \leq %	CV \leq %
1 - 50	50	1.8	0.8
	25	2.5	1.2
	5	8.0	5.0
10 - 200	200	1.0	0.3
	100	1.5	0.4
	20	4.0	1.5
40 - 1000	1000	1.0	0.2
	500	1.5	0.3
	100	3.0	1.0

Table 4: Accuracy Table 1-channel Liquid End

8-channel Liquid End

Volume range [μl]	Partial volume [μl]	A \leq %	CV \leq %
1 - 50	50	1.5	1.0
	25	2.0	1.6
	5	8.0	20.0
20 - 300	300	1.2	0.4
	150	1.6	0.6
	30	5.0	2.5
40 - 1000	1000	1.2	0.3
	500	1.6	0.5
	100	4	2

Table 5: Accuracy Table 8-channel Liquid End

1.6 Commissioning

1.6.1 Installation

 **DANGER**

Explosion hazard

Never operate the device in potentially explosive atmospheres.

The surface must be level and have a minimum carrying capacity of 25 kg for the Liquid Handling Station and 35 kg for the Liquid Handling Station flow.

Verify that the device is installed horizontally level with a spirit balance.

The recommended minimum installation space is H 700 mm / 750 mm x W 800 mm x D 600 mm. These dimensions guarantee a distance of approx. 100 mm of all sides of the device to the wall, thereby allowing unobstructed air circulation and protecting the device from overheating. The dimensions also allow 160 mm of free space on top, which is required to open the front door.

Select a location that meets the following operating conditions:

- + 15 °C to + 35°C
- max. 95 % relative humidity at 30 °C

Moreover, provide for sufficient space for a monitor and a PC or for a laptop.

NOTICE! Keep temperature differences and changes in humidity levels to a minimum during setup and transport of the machine. Condensation can lead to malfunctions and defects in the Liquid Handling Station or the computer connected to it.

Prior to commissioning, the Liquid Handling Station and the computer must first adapt to the ambient temperature, which may require several hours under certain conditions.

NOTICE

Risk of material damage

In order to prevent material damage, protect the device from:

- Dust and drafts.
 - Shocks and vibrations.
 - Electromagnetic fields (e.g. motors).
 - Aggressive vapors (e.g. corrosive media).
 - Dripping and splashing water. Do not place any vessels containing liquids on or immediately next to the device.
 - Exposure to direct sunlight or radiation from radiators.
 - Moreover, ensure that air can circulate freely around the device.
-

1.6.2 Connection

⚠ DANGER

Danger to life from electric current

Touching live parts poses a danger to life due to electric shock

- Touching live parts poses a danger to life due to electric shock
- Ensure at all times that the housing and covers are intact and closed.
- Only connect the device to sockets with a protective contact.
- Before connecting the device, check whether the mains voltage agrees with the admissible operating voltage. An incorrect mains voltage may result in the destruction of the device.
- Before switching on the device, ensure that the power cable and connections are intact.
- The device may only be opened by qualified persons authorized by the manufacturer.
- Before opening or removing covers, switch off the device and pull out the power plug and the USB plug.
- When using the *Liquid Handling Station flow*, the FlowBox must also be switched off and the DC 24V plug and the USB plug pulled out.
- Keep humidity away from live parts in order to prevent short circuits.

NOTICE

User qualifications

In order to ensure that the device functions properly and no hazards are posed to the user, you must have the device installed and initial commissioning carried out by the manufacturer. Moreover, all users of the device must have received training from the manufacturer or an authorized dealer prior to commissioning.

NOTICE

USB interface

The connection is via USB cable; communication between the PC and the Liquid Handling Station and Flowbox takes place via a virtual COM interface. For this reason, both the Liquid Handling Station and the Flowbox appear in the Windows Device Manager under Ports | Communication Port (COM).

The manufacturer must be contacted prior to the decommissioning, transport, and disposal of the device.

1.6.3 Connecting the Liquid Handling Station

- a. Connect the *Liquid Handling Station* to the PC with the USB cable.
- b. Connect the Liquid Handling Station to the power supply with the included power cable.
- c. Switch on the power button on the back.

1.6.4 Connecting the Liquid Handling Station flow


- a. Connect the *Liquid Handling Station* to the PC with the USB cable.
 - b. Connect the Liquid Handling Station to the power supply with the included power cable.
 - c. Switch on the power button on the back.
 - d. Connect the *FlowBox* to the PC with the USB cable.
 - e. Connect the *FlowBox* with the included power adapter.
- ➔ The LED on the back indicates that the power supply is on.

1.6.5 Connecting the shaking and temperature control module Teleshake 95

NOTICE

Component damage due to missing potential equalization

When installing the Teleshake 95 components, bear in mind that they may be damaged by electrostatic discharge. In order to equalize their electrical potential, touch a metallic surface (e.g. radiator). Proceed with the installation only after you have done this. Also observe the INHECO instruction manual: <https://www.inheco.com/service/downloads.html>.

Ill. 5: Teleshake 95 installation

- a. Pull the power plug of the *Liquid Handling Station* and open the front door.
- b. Remove the blind plug at the opening on the back of the device.
- c. Set the Teleshake 95 to one of the positions P2 ... P4 and run the Sub-D connector through the opening. **NOTICE! The positions P2 ... P4 are suitable due to their proximity to the opening on the back of the device.**
- d. Connect the Sub-D connector to a slot module. **NOTICE! The color of the shrink sleeve on the connector must match the color coding next to the slot module on the control unit.**
- e. Connect the USB cable of the *Liquid Handling Station*, *FlowBox* (if available) and *MTC* or *STC* control unit to the USB hub and connect the USB hub to your laptop/PC.
- f. Configure the Teleshake 95. Also observe the INHECO instruction manual: <https://www.inheco.com/service/downloads.html>.

1.7 Operation

1.7.1 Liquid Handling Station

NOTICE

Equip the work table correctly

- Make absolutely sure that the labware is placed on the work table as entered in the software/method.
- In order to compensate for height differences, it is absolutely necessary to use the adapters recommended by the software.
- Non-observance poses the risk that the Liquid End will strike the labware and become damaged.
- Only use height adapters in the designated direction. Do not apply force.

- a. Switch on the device with the on/off switch on the front right of the device.
- b. Start the software on the PC. A connection is established between the device and the PC.
- c. During initial commissioning or after creating a new workbench, select the device type that you acquired under System | Settings | Device Manager | Device Type.
- d. Create a method or load an already existing method. This procedure is described in detail in the separate user manual for the Liquid Handling Station operating software.
- e. Open the front door and equip the work table with the labware that is required for the method and that was selected in the software.
- f. Close the front door and start the program.

⚠ CAUTION

Risk of crushing

Only reach into the work area when the device is at a standstill!

NOTICE**Monitored front door**

In delivery condition, the program can only be started when the front door is closed. If the front door is opened during a program, the operation is interrupted once the current work step has been completed. Once the front door is closed again, the program continues from the point of interruption.

1.7.2 Liquid Handling Station flow**NOTICE****Unusual operating noise**

Pay attention to changing operating noises while operating the FlowBox. If the operating noise increases with time, this may indicate damage. Contact BRAND.

NOTICE**Reducing turbulence**

In order to reduce turbulence, the front door should always be opened slowly and as rarely as possible and the protective screen inserted.

NOTICE**Reporting FlowBox**


The FlowBox can be switched on both manually and by means of the software. If the FlowBox activity should appear in the report, the FlowBox must be switched on by means of the software!

The FlowBox can only be switched on once the Liquid Handling Station flow is ready for operation and the FlowBox has been connected to the power supply with the included power adapter.

- a. Decontaminate the interior of the Liquid Handling Station flow according to the specifications provided in the chapter "Maintenance and cleaning."
- b. Decontaminate the labware used, the adapters, the Liquid Ends, and the protective screen.
- c. Insert the protective screen.
- d. Close the front door.
- e. Switch on the Liquid Handling Station flow using the on/off switch.
- f. Switch on the FlowBox using the on/off switch.
- g. Allow the FlowBox to run for 5 minutes in order to flush/filter all the air in the work area of the pipetting robot.
- h. Start your method using the software.
- i. Equip the work table according to the assignment request of the software.
- j. Close the front door.

k. After completing the method and removing your samples, switch off the FlowBox.

1.7.3 Liquid End

Ill. 6: Exchanging a Liquid End (example shows a 1-channel Liquid End)

NOTICE

1-channel or 8-channel Liquid Endl

The figure illustrates a 1-channel Liquid End. An 8-channel Liquid End is changed using the same procedure.

NOTICE

Exchanging a Liquid End while the Liquid Handling Station is off

If you want to change the *Liquid End* while the Liquid Handling Station is switched off, it has to be pulled down while pressing the decoupling button.

NOTICE

Different Liquid End in method sequence

If you require different Liquid Ends for the execution of a method, a message on the screen will prompt you to change the Liquid End.

1.7.4 Shaking and temperature control module Teleshake 95

The INHECO *Teleshake 95* is a shaking and heating module and is placed on one of the positions P2 ... P4 of the *Work Table*. The connection cable is guided through an opening in the rear wall and connected to the INHECO STC or MTC control unit via color-coded connectors. After completing the installation, the *Teleshake 95* must be activated in the controller and set up for the respective method in the *Device Commands*.

See:

- To install the *Teleshake 95*, see Installation.
- To activate the *Teleshake 95*, see Device Manager.
- To set up the *Teleshake 95*, see Shake Command or Heat Command.

1.8 Maintenance/Cleaning

1.8.1 Device

Immediately wipe off any splashes, drops or large quantities of spilled liquids with an absorbent cloth.

Remove any contaminants (e.g. dirt, dust) with a soft, clean cloth.

If necessary, use a customary, neutral cleaning agent.

The Liquid Ends are autoclavable for 20 minutes at 121°C.

The adapters are not autoclavable! They can be cleaned where necessary with ethanol or disinfectants.

Decontamination

⚠ CAUTION

Switching off the device and disconnecting the power

Before carrying out cleaning, decontamination or maintenance, switch off the Liquid Handling Station and disconnect it from the power supply. When using the Liquid Handling Station flow, the FlowBox must also be switched off and the DC 24V plug pulled out.

Wipe off the device and work area with a lint-free cloth.

To disinfect or decontaminate the device, we recommend using 70% (v/v) ethanol, 3% ... 4% sodium hypochloride or alcohol-based surface disinfectant cleaner, e.g. Pursept® -A Xpress disinfectant spray.

NOTICE

Damage to device during cleaning and decontamination

When spraying, ensure that

- the motor control unit does not come into contact with cleaning agents and other liquids.
- no cleaning agents enter the guide slot.
- no cleaning agents get on the discharge filter and the intake filter when using the Liquid Handling Station flow.
- If liquid enters the device, immediately disconnect the device from the power supply and contact BRAND Service.

A UV decontamination of the Liquid Ends and adapters is possible. If another method or cleaning agent is used besides those mentioned here, please clarify with the manufacturer whether it can potentially damage the device.

▲ WARNING

Cleaning in medical laboratories

Special care must be taken when handling infectious material (see Danger from harmful substances or infectious media, p. 60). Observe the relevant safety regulations.

Shipment for service or repair

If the device has to be shipped to the manufacturer, a decontamination must be carried out and documented beforehand. For safety reasons, only decontaminated devices can be serviced and repaired.

For this purpose, fill out the “Declaration of Health Safety” and send it together with the device to the manufacturer. Pre-printed forms can be requested from the manufacturer or are available for download at www.brand.de.

Before sending the device, request the special packaging from the manufacturer and pack the Liquid Handling Station properly in order to prevent damage during transport.


NOTICE

Transporting hazardous materials without approval is prohibited by law.

1.8.2 Liquid Ends

1.8.2.1 1-channel Liquid End

In order to ensure their proper functioning, the Liquid Ends should be serviced in regular intervals and cleaned where necessary.

Ill. 7: 1-Channel Liquid End

- | | |
|------------------------------|------------------------------|
| 1 Piston unit | 4 Lower part of ejector unit |
| 2 Upper part of ejector unit | 5 Shaft |
| 3 Seal | 6 Liquid End |

Maintenance

- a. Decouple the Liquid End: Press the decoupling button with one hand and simultaneously remove the Liquid End with the other hand.
- b. Check the tip mounting cone at the tip of the pipette shaft (5) for damage.
- c. Inspect the piston at the lower end of the piston unit (1) and the seal (3) for damage or contamination. The liquid end must be dismantled for this purpose. Replace the piston if necessary.
- d. Check the tightness of the liquid end.

NOTICE

For the check, we recommend the BRAND PLT leak testing unit. If the Liquid End is not leak-tight, please contact the Support department of the manufacturer.

Disassembly and cleaning

- a. Decouple the liquid end from the motor control unit.
- b. Unscrew the retention bushing (1). During this, the piston unit remains connected to the retention bushing.
- c. Remove the spring with seal (3).
- d. Unscrew the upper part of the ejector unit (2) from the pipette shaft.
- e. Pull out the shaft (4) from the lower part of the ejector unit (5).
- f. Clean the illustrated parts with soap solution or isopropanol. Afterwards, rinse with aqua dest.
- g. Allow the parts to dry or dry at max. 120 °C.
- h. Reoil pistons with a very thin coat of the supplied silicone oil.
- i. Reassemble the cooled parts in reverse order.
- j. Only hand-tighten the retention bushing (1) and the upper part of the ejector unit (2).

NOTICE

All of the individual components shown in the figure can be procured as spare parts. Use only original spare parts!

1.8.2.2 8-channel Liquid End

In order to ensure their proper functioning, the Liquid Ends should be serviced in regular intervals and cleaned where necessary.

Maintenance

- a. Decouple the Liquid End.
- b. Check shafts, pistons, and seals for damage and contamination.
- c. Check that the Liquid End is leak-tight.

NOTICE

Check the circulation

For the check, we recommend the BRAND PLT leak testing unit. Use only original spare parts!


Disassembly and cleaning

- Individual shafts, pistons, and shaft/piston support bars (these parts only) may be cleaned with soap solution or isopropanol. Afterwards, rinse with aqua dest.

- Allow the parts to dry completely or dry at max. 120 °C and allow to cool. Liquid residues in the shafts result in deviations in accuracy.
- Reoil pistons with a very thin coat of the supplied silicone oil. For the central guide axis, only use the prescribed fluorine static grease!


1.8.2.2.1 Dismantle piston unit (up to 300 µl MC8)

Check that the Liquid End is leak-tight.


- Decouple the Liquid End
- Rotate the two locks by 90° and remove the lower part of the housing.

Unscrew the shaft

- Apply mounting tool on individual shaft.
- Unscrew the shaft.


Remove seal

- Push the piston unit all the way down.
- Remove and inspect the seal and clean or change where necessary.

After removing the shaft, the seal is located either in the shaft or on the piston.

Mount the seal

- a. If necessary, reoil the piston slightly with the supplied silicone oil. Only use the recommended silicone oil.
- b. Slide the seal onto the piston with the flat side upwards.
- c. Tighten the cleaned or new shaft using the supplied mounting tool.

Couple the Liquid End


- a. Reassemble the Liquid End, mount the lower part of the ejector unit, and check for leak-tightness.
- b. After checking the Liquid End, couple it onto the motor control unit again.

NOTICE

For the check, we recommend the BRAND PLT leak testing unit. If the Liquid End is not leak-tight, please contact the Support department of the manufacturer.


1.8.2.2 Dismantle piston unit (1000 µl MC8)

Remove housing of the liquid end


- a. Decouple the Liquid End.
- b. Pull out sliding clamps of the housing cover laterally and remove the housing.

Remove shaft fastening


- a. Pull out the shaft fastening laterally.

Disassemble the shaft and piston

- a. Carefully press the shaft out forward. Take the piston directly with you and slowly pull it out of the holder.
- b. To re-grease, pull the piston out of the shaft. Only all piston-shaft units of a Liquid End the can be replaced.

Regrease piston seal


NOTICE

Apply silicone grease only very thinly to the seal with a brush.


- a. Pull the piston with seal out of the shaft.
- b. Dip the brush into grease and wipe well along the wall of the vessel.
- c. Hold brush against seal and turn piston with seal 1-2 turns.
- d. Insert the piston with seal back into the shaft.

Assemble shaft and piston

- a. Push the shaft into the holder and carefully press the piston back into its original position (if necessary, press the piston spring slightly upwards). The shafts must then be aligned again.
- b. Continue with all piston-shaft units.

Assembling the pipetting unit.


- a. Slide in the shaft fastening.
- b. Slide housing onto pipetting unit.
- c. Push the slide fasteners back in.

1.8.2.2.3 Make piston unit freely movable (1000 µl MC8)


After long storage time it is possible that the piston unit of the 1000 µl MC8 cannot be moved by the motor control unit. In this case proceed as follows: Proceed as follows:

Remove housing of the liquid end

- a. Decouple the Liquid End.
- b. Pull out sliding clamps of the housing cover laterally and remove the housing.

Free movement of piston unit

- a. Move piston unit several times by hand up to the upper and lower stop.
 - b. Slide housing onto pipetting unit.
 - c. Push the slide fasteners back in.
- ↳ The 1000 µl MC8 is again free to move.

1.9 Ordering Information

1.9.1 Liquid Handling Station and accessories

1.9.1.1 Software

Description	Order No.
Single user license	7094 90
Upgrade to new version	7094 92

1.9.1.2 Pipetting robot

1.9.1.2.1 BRAND Liquid Handling Station

including motor control unit, operating software, operating manual, mains cable, USB cable, documentation, and on-site training.

Order No. 7094 00

1.9.1.2.2 BRAND Liquid Handling Station flow

including motor control unit, operating software, operating manual, 2 mains cables, 2 USB cables, 24 V DC universal power adapter for FlowBox, protective screen, documentation, and on-site training.

Order No. 7094 02

1.9.1.3 Teleshake 95 shaking and heating module

Description	Order No.
Teleshake 95 shaking and heating module with permanently mounted centering plate for the Liquid Handling Station	7094 70

1.9.1.3.1 Controller for Teleshake 95

1.9.1.3.1.1 Single TEC Control

Designation	Order No.
Single TEC Control + 1 black slot module	7094 72

1.9.1.3.1.2 Multi TEC Control

Description	Order No.
Multi TEC Control + 2 black slot modules	7094 74

1.9.1.3.2 Adapter for Teleshake 95

Designation	Order No.
Universal adapter for flat bottom plates. This universal adapter is included with the Teleshake 95.	7094 76
Universal adapter for 96-well PCR plates, low profile	7094 77
Universal adapter for 96-well PCR plates	7094 78
Universal adapter for 384-well PCR plates	7094 79
Adapter for 96-well BRAND F- and C-bottom plates	7094 80
Adapter for 96-well BRAND V-bottom plates	7094 81
Adapter for 96-well BRAND U-bottom plates	7094 82

1.9.1.4 USB hub

Description	Packaging unit	Order No.
USB Hub connection extension	1	7094 89

1.9.1.5 Liquid Ends and Accessories

1.9.1.5.1 1-channel Liquid End

Volume	Description	Order No.
1 - 50 µl	SC-50	7094 10
10 - 200 µl	SC-200	7094 13
40 - 1000 µl	SC-1000	7094 16

1.9.1.5.2 8-channel Liquid End

Volume	Description	Order No.
1 - 50 µl	MC-50	7094 20
20 - 300 µl	MC-300	7094 23

Volume	Description	Order No.
40 – 1000 µl	MC-1000	7094 26

1.9.1.5.3 Liquid End holder

For use with 1-channel and 8-channel Liquid Ends, pkg unit 1 pc.

Description	Order No.
for 3 Liquid Ends	7094 63
for 5 Liquid Ends	7094 65

1.9.1.5.4 Grease and Lubricants

Designation	Order No.
Silicone oil, for pistons, 1-channel Liquid Ends	705502
Silicone grease, for pistons, 8-channel Liquid Ends	703677
Fluorine static grease, for guide axis, 8-channel Liquid Ends	703678
PFPE Piston grease, for 1000µl 8-channel Liquid Ends	61613

1.9.1.6 Labware

1.9.1.6.1 Labware-Adapter and Racks

Pkg unit 1 pc.

Description	Material	Order No.
Height adapter, 60 mm	Aluminum	7094 30
Height adapter, 30 mm	Aluminum	7094 32
Tip adapter	Aluminum	7094 34
Cuvette adapter	POM	7094 36
Reservoir rack for four 40 ml reservoirs (701460, 701462)	Aluminum	7094 43
PCR adapter 96-well	POM	7094 46
PCR adapter 384-well	POM	7094 48
PCR adapter 384-well, optimized for BRAND plates (781345, 781347 and 781358)	POM	7094 49
Micro-tube rack 1.5 ml	POM	7094 50
Micro-tube rack 0.5 ml	POM	7094 52
Micro-tube rack 5 ml	POM	7094 53
12 x 75 mm tube rack	POM	7094 55

1.9.1.6.2 PCR cooling rack

Description	Order No.
PP. Indirect cooling, for 0.2 ml single tubes, strips of 8 and 12, and 96-well PCR plates. Changes from violet to pink at 7°C. Pkg unit 2 pc	7094 56

Description	Order No.
Pkg unit 2 pc	

1.9.1.6.3 Waste box

Description	Order No.
including reservoir for residual liquids, pkg unit 5 pc.	7094 58

1.9.1.6.4 Reagent reservoirs


PP, autoklavierbar.

Volume	Design	Pkg unit	Order No.
220 ml	pyramid bottom, 96-well	10	7014 50
12 x 6 ml	pyramid bottom, low profile	10	7014 52
4 x 60 ml	pyramid bottom	10	7014 54
6 x 40 ml	pyramid bottom	10	7014 56
50 ml	flat bottom, low profile	10	7014 58
40 ml	v-shaped bottom, low dead volume	24	7014 60
40 ml	>v-shaped bottom, low dead volume	24	7014 62


1.9.1.6.5 Passive cooler Cooling Block

The Cooling Block can also be used with suitable disposables from other manufacturers. However, the suitability test must be carried out by the user themselves. The “Cooling Block” passive cooler is patent pending.

Cooling Block for microcentrifuge tubes, 0.5 ml

Order No. 7095 10	Suitable BRAND labware
	<p>BRAND Tube 0,5 ml: 780750, 780755, 780760, 781310, 781311, 781312, 781313, 781314</p> <p>BRAND PCR single tubes, 0.5 ml: 781310, 781311, 781312, 781313, 781314</p> <p>BRAND microcentrifuge tubes 0,5 ml: 780700, 780710, 780730, 780750, 781313, 781314</p>


Cooling Block for microcentrifuge tubes, 1.5 ml

Order No. 7095 11	Suitable BRAND labware
 <p>The image shows a rectangular cooling block with a white base and a blue top. The top surface is a grid of 24 circular wells, arranged in 4 rows and 6 columns. The block is shown from a three-quarter perspective.</p>	<p>BRAND microcentrifuge tubes 1.5 ml: 780500, 780502, 780400, 780521, 780522, 780523, 780524, 780525, 780540, 780505, 780701, 780702, 780711, 780712, 780751, 780752, 780731, 780732, 780756, 780757, 780761, 780762</p> <p>BRAND microcentrifuge tubes 2.0 ml: 780546, 780550, 780703, 780704, 780713, 780714, 780733, 780734, 780753, 780754, 780758, 780759, 780763, 780764</p>

Cooling Block for microcentrifuge tubes, 5.0 ml


Order No. 7095 12	Suitable BRAND labware
 <p>The image shows a rectangular cooling block with a white base and a blue top. The top surface is a grid of 12 circular wells, arranged in 3 rows and 4 columns. The block is shown from a three-quarter perspective.</p>	<p>No available BRAND microcentrifuge tubes</p>

Cooling Block for 96 well PCR plates

Order No. 7095 13	Suitable BRAND labware
 <p>The image shows a rectangular cooling block with a white base and a blue top. The top surface is a grid of 96 small circular wells, arranged in 8 rows and 12 columns. The block is shown from a three-quarter perspective.</p>	<p>BRAND PCR 96 plates: 781350, 781354, 781357, 781364, 781365, 781366, 781367, 781368, 781369, 781371, 781372, 781373, 781374, 781375, 781376, 781377, 781378, 781400</p> <p>BRAND 24-well PCR plates: 781411, 781412</p> <p>BRAND 48-well PCR plates: 781415, 781416</p>

Order No. 7095 13	Suitable BRAND labware
	<p>BRAND 8-tube PCR strips, 0.2 ml:: 780500, 780502, 781315, 781316, 781317, 781320, 781321, 781322, 781323, 781324, 781325, 781326, 781327, 781330, 781332, 781333, 781340, 781341, 781342, 781343, 781344, 781334</p> <p>BRAND 12-tube PCR strips, 0.2 ml: 781280, 781284, 781290</p> <p>BRAND PCR single tubes, 0.2 ml: 781300, 781301, 781302, 781303, 781304, 781305</p>

Cooling Block for 384 well PCR plates

Order No. 7095 14	Suitable BRAND labware
	<p>BRAND 384 well PCR plates: 781345, 781347, 781348, 781358</p>

1.9.1.7 Tips

Selection

Use the following compatibility table to find the right tip for your Liquid Ends:

Liquid Ends		Robot tips [µl]				Robot filter tips [µl]			
Order No.	Des. in software	1- 50	10- 200	10- 300	40- 1000	1- 20	10- 100	10- 200	40- 1000
709410	1 channel 50 µl SC-50 µl	✓	—	—	—	✓	—	—	—
709413	1 channel 10 µl SC-200 µl	—	✓	✓	—	—	✓	✓	—
709416	1 channel 40 µl SC-1000 µl	—	—	—	✓	—	—	—	✓
709420	8 channel 1-50 µl MC8-50 µl	✓	—	—	—	✓	—	—	—

Liquid Ends		Robot tips [μ l]				Robot filter tips [μ l]			
709423	8 channel 20-300 μ l MC8-300 μ l	—	✓	✓	—	—	✓	✓	—
709426	8 channel 40-1000 μ l MC8-1000 μ l	—	—	—	✓	—	—	—	✓

Robot tips

Volume	Packing unit	Order No. non-sterile	Order No. sterile
1-50 μ l	TipRacks à 96	733146	733166
10-200 μ l	TipRacks à 96	734148	734168
10-300 μ l	TipRacks à 96	734150	734170
40-10000 μ l	TipRacks à 96	734152	734172

Robot filter tips

Volume	Packing unit	Order No. non-sterile	Order No. sterile
1-20 μ l	TipRacks à 96	733646	733666
10-100 μ l	TipRacks à 96	734650	734670
10-200 μ l	TipRacks à 96	734652	734672
40-10000 μ l	TipRacks à 96	734654	734674

1.9.1.8 Leak detector for liquid ends.


Description	Order No.
PLT unit, pipette leak detector	7039 70

1.9.2 Spare parts

1.9.2.1 Liquid End 1-Kanal 200 -1000 μ l

Packing unit 1 piece

The appearance and dimensions of the spare parts correspond to the respective nominal volume (fig. 200 μ l Liquid End)

Volume	A	B*	C	D	E
1 - 50 µl	7096 02	-	7096 08	7096 12	7096 24

*not a standard spare part, please contact the manufacturer.

Piston and seal sets

Volume	Order No.
10 - 200 µl	7096 04
40 - 1000 µl	7096 06

1.9.2.2 Liquid End 8-Channel 50 - 300 µl

Packing unit (A) 8 pieces, PU (B) 3 pieces

The appearance and dimensions of the spare parts correspond to the respective nominal volume.

Ill. 8: Shaft and seal

- 1 Seal
- 2 Shaft

Volume	A*	B
1 - 50 µl	7096 30	7033 43
20-300 µl	7096 32	7033 46

*incl. 8 gaskets and 1 assembly wrench


Liquid End 20-300µl additionally with pressure ring.

NOTICE

Shaft change

If a socket change becomes necessary, all sockets of a multi-channel liquid end must always be replaced!

1.9.2.3 Liquid End 8-Channel 1000 µl

1 Piston with seal

2 Shaft

Volume	Designation	Order No.
40 – 1000 µl	8 x Flask unit for 1000 µl MC8 Liquid End	709634

2 Glossary

A

Adapter

Adapters are used to compensate for height differences between the different labware types. They can be used to achieve shorter paths and a faster implementation of the method. The distinction is made between height adapters, tip adapters (for tips), micro-tubes, and PCR racks. Overview of all available adapters, see Ordering Information, p. 84.

C

Command

Refers to a single sequential step in the program sequence of a method. In the operating software, four different commands are generally available: Transfer, Wait, Mix, and Break Command. Additional device commands are available depending on the device type and accessories.

Cooler

Refers to a device used to cool samples.

Content structure

The content structure of a report comprises all the elements between the report header and footer, i.e. the actual (method-specific) content of the report.

D

Destination

Labware or target vessel into which liquid from another vessel (source) is transferred.

Destination well

Single well of labware defined as a destination.

F

FlowBox

Permanently installed filter device that purifies the air in the work area of the Liquid Handling Station flow.

Filling volume

Maximum volume of a labware or an individual well available for filling liquid.

G

Graphical view

Area in software that graphically displays the work table (with selected labware, filling, etc.).

H

Heater and shaker

The term *Heater* and *Shaker* refers to the device *Teleshake 95*, which can be used to heat and/or shake samples over an adjustable period of time. The device is placed at a position in the *Work Table* and controlled via the *Heat/Shake* command. The device can also be controlled outside of the *Liquid Handling Station Software* in order, for example, to heat samples over a period of time longer than the method sequence.

L

Labware

Collective term for the microplates, PCR plates, deep-well plates, inserts, tube racks, tubes, tips, reservoirs, and waste box that can be used with the Liquid Handling Station.

Liquid End

1- or 8-channel pipetting unit coupled to the motor control unit in order to allow the execution of liquid transfers with the Liquid Handling Station.

Liquid type

Definition of various liquid parameters, such as aspiration and dispensing speed.

M

Master layout

The master layout of a report comprises the report header and footer, as well as general settings for the page layout, e.g. page margins, font size and style in the report. **Method**

Procedure for the loading of the work table, the definition of commands, and the execution of the method. You can find a more detailed explanation of the term “method” in the chapter “Method Structure”.

Minimum immersion depth

(= minimum tip immersion depth). Distance of the lower end of the tip to the liquid surface.

This distance guarantees correct liquid aspiration and prevents, for example, accidental air intake.

Minimum bottom distance

(= minimum bottom distance). Refers to the shortest distance between the vessel bottom and the lower end of the tip moved into the labware during immersion. This distance is used to prevent crashes and to ensure a correct liquid aspiration.

Multidispense

Pipetting mode in which a volume aspirated from a source is dispensed into a destination in partial steps.

P

Pipette

Pipetting mode in which a previously defined volume is aspirated from a source and then completely dispensed immediately thereafter into a destination.

Pipetting step

The aspiration of a previously defined volume from a source and dispensing of this volume (as a whole or in partial steps).

Position

Individual recess on the work table of the Liquid Handling Station on which labware can be positioned. The work table of the Liquid Handling Station has eight positions. Seven of these can be equipped as desired, while one position is designated for the waste box.

Prewet

Designation for the preliminary conditioning of the tips through liquid aspiration from the source well and subsequent liquid dispensing into the source well.

Priority

Setting option for liquid names in the well content. If the liquid name is given a high priority, this is retained for all methods, i.e. the transfer history of a certain liquid can be tracked.

R**Remaining volume**

(= remaining volume). Volume that must remain in the vessel to allow a defined and reliable intake. It depends on the minimum immersion depth and the minimum bottom distance.

Restart behavior

Setting option for liquid volumes. If the method is restarted, the volume remaining from the feed line is used as the starting volume.

S**Sample selection**

Setting in the Execute tab that enables the reduction of the existing methods to the actual number of accumulating samples.

Settings

Settings that can be made for the labware at the individual positions and for the different commands.

Shaker

see Heater/Shaker

SMTP settings

Settings for the login of the e-mail service.

Source

Labware or source vessel from which liquid is removed in order to be transferred to another vessel (destination).

Source well

Individual well of labware defined as a source.

T**Tab**

Refers to individual tabs on a menu tab. Example: Method tab in the main menu tab.

Transfer

A single transfer includes all the pipetting steps executed between one or more sources and one or more destinations in a (Transfer) command.

Transfer volume

Liquid volume that is transferred from the source to the destination or into the individual wells of the destination.

V**Distribution pattern**

Specification of the source wells and destination wells.

W

Well content

In general: The content of one or more wells. In the operating software: Settings (e.g. liquid volume, name, liquid type) that apply to the content of a fillable labware.

Workbench

All relevant data required for using the Liquid Handling Station Software is found in the Workbench. (E.g. labware definitions, methods, user settings, etc.).

Work table

Work area of the Liquid Handling Station whose eight positions can be virtually loaded using the software. The definition of the commands is followed by the physical loading of the work table.

Z

Z-tracking

Tracking of the liquid end on the Z-axis during liquid aspiration and dispensing. The Z-tracking corresponds to a downward movement during liquid aspiration and an upward movement during liquid dispensing. The starting point of aspirating and dispensing can be user-defined if necessary.

997496 | Printed in Germany | 0522