[Gebrauchsanleitung](#) | [Operating manual](#) | [Mode d'emploi](#) | [Instrucciones de manejo](#) | [Istruzione](#) | [Instruções de utilização](#) | [操作手册](#) | [Руководство по эксплуатации](#) | [사용 지침](#) | [Bruksanvisning](#) | [使用説明書](#) | [Használati utasítás](#) | [Návod k použití](#) | [Gebruiksaanwijzing](#) | [Instrukcja użytkowania](#) | [Kullanım Talimatları](#)

Transferpette® pro

Pipettierhelfer | Pipette controllers

Impressum

BRAND GMBH + CO KG

Otto-Schott-Str. 25
97877 Wertheim (Germany)

T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

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Table of contents

1	Using the operating manual	51
1.1	Signal words and their meaning	51
1.2	Presentation of descriptions of actions	51
2	Safety regulations	52
2.1	General safety regulations	52
2.2	Target group	53
2.3	Purpose	53
2.4	Use	54
2.5	Improper use	54
2.6	Foreseeable misuse	54
2.7	Operating Limitations	54
2.8	Limitations of Use	54
2.9	Operating exclusions	55
3	Scope of delivery	56
4	Functional and operational components	57
5	Pipetting	61
6	Storage	65
7	Checking the volume	66
8	Accuracy table	68
9	Adjusting – User and factory adjustment	70
9.1	Permanent factory adjustment: Easy Calibration	70
9.2	Temporary adjustment: User Adjustment	72
10	Disinfection/Autoclaving	76
10.1	Autoclaving	76
10.2	UV disinfection	76
10.3	PE filter	76

11 Maintenance	78
11.1 Cleaning – volumes up to 1,000 μ L.....	79
11.2 Cleaning – volumes of 2.5, 5, and 10 ml.....	80
12 Troubleshooting	81
13 Product markings	82
14 Order Information	83
15 Repairs	87
15.1 Sending for repair	87
16 Calibration service	89
17 Information about your laboratory device	90
18 Warranty	91
19 Disposal	92

1 Using the operating manual

- Please carefully read the operating manual before using the device for the first time.
- Keep the operating manual in an easily accessible place. It is part of the instrument.
- Be sure to include the operating manual if you transfer possession of this device to a third party.

1.1 Signal words and their meaning

⚠ WARNING

or **⚠WARNING! ...**

The signal word WARNING indicates a hazardous situation that, if not avoided, could result in death or serious injury.

⚠ CAUTION

or **⚠CAUTION! ...**

The signal word CAUTION indicates a hazardous situation that, if not avoided, could result in minor or moderate injury.

NOTICE

or **NOTICE ...**

The signal word NOTICE is used to address practices not related to physical injury. Example: Possible property damage.

1.2 Presentation of descriptions of actions

Format	Meaning
1. Task	Indicates a task.
a., b., c.	Indicates the individual steps of a task.
>	Indicates a prerequisite for a task.
⇒	Indicates a result of a completed task.

2 Safety regulations

2.1 General safety regulations

Please read carefully!

The Transferpette® S can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

- Every user must read and observe the operating manual enclosed with this laboratory device before and during use. The device may only be used by trained and qualified personnel.
- Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves.
- When working with infectious or hazardous samples/media (e.g. hazardous substances), the general safety rules in the laboratory must be observed and regulations for handling the samples/media must be followed. The information provided by the media manufacturer (e.g. safety data sheets) must be observed.
- Use the laboratory device only for pipetting or dispensing media within the defined limitations and restrictions of use. Comply with the operating exclusions.
- If working with inflammable media, make sure to take the necessary precautions to avoid the buildup of static electricity (e.g., do not dispense into plastic vessels and do not wipe instruments with a dry cloth). Do not use the device in potentially explosive atmospheres. If in doubt, contact the manufacturer or supplier.
- Always check that the instrument is in proper working condition before use. If device malfunctions are indicated (e.g. sluggish pistons, leaks or power supply faults), stop working with the device immediately and refer to the "Troubleshooting" section in the operating manual. Contact the manufacturer, if necessary.

- Always perform work in a manner that does not endanger users or other people. Avoid splashes. Only use suitable vessels. Never use unnecessary force when operating, cleaning or servicing the laboratory device.
- If the laboratory device is supplied with power through a power adapter, batteries or accumulators, the components and connections on the device must be regularly checked for proper working condition. Do not operate the laboratory device and its accessories in an unprotected, humid or wet environment.
- Do not make any technical modifications. Use only original replacement parts of the manufacturer; even power adapters or accumulators of identical size and specification from other manufacturers are not permitted. Do not disassemble the laboratory device and its accessories (e.g., power adapters, cables, racks, accumulators or batteries) any further than described in the user manual!
- Only autoclave the laboratory device when this is permissible according to the user manual.

2.2 Target group

The operating manual is intended for users who use the laboratory instrument in the course of their professional activities. Users are familiar with the typical safety regulations and working methods in laboratories and have been trained accordingly. They can recognize potential hazards and protect themselves from them. The operating manual assumes this expertise and does not replace basic laboratory training or specific safety training.

2.3 Purpose

This is an air displacement pipette for pipetting liquids of medium density and low to medium viscosity.

2.4 Use

Use the laboratory device Transferpette® S only for pipetting or dispensing liquids within the defined limits of use.

2.5 Improper use

Various risks may arise if the laboratory instrument is used improperly. These risks include: inaccurate liquid delivery, damage to the laboratory instrument, and the risk of contamination, infection, and injury from contact with the pipetted media.

Any use other than for pipetting or dispensing liquids within the defined operating limits is considered improper use.

2.6 Foreseeable misuse

A typical misuse is pipetting or dispensing liquids with too high viscosity or using unsuitable tips.

2.7 Operating Limitations

Volume accuracy may be affected by liquid density, special tip shapes, or deviations from ambient temperature.

However, the temporary user adjustment enables you to correct any resulting deviations and improves accuracy under conditions that differ from the factory setting (aqueous medium, ISO 8655). See Temporary adjustment: User Adjustment, p. 72 .

2.8 Limitations of Use

The pipette is used for dispensing liquids within the following limits:

- Operating temperature from +15 °C to +40 °C (59 °F to 104 °F).
Additional temperatures upon request
- vapor pressure up to 500 mbar
- Viscosity: 260 mPa s

With viscous liquids, the dispensing speed may need to be adjusted.

2.9 Operating exclusions

The user must verify the suitability of the instrument for the intended purpose, as aggressive liquids and their vapors can damage the instrument (corrosion!). The instrument cannot be used for the following liquids:

- for liquids with very high steam pressure
- Liquids that corrode the following materials:
 - Fluoroelastomer rubber (FKM)
 - Polyamide (PA)
 - Polycarbonate (inspection window)
 - Polyetheretherketone (PEEK)
 - Polyphenylene sulfide (PPS)
 - Polypropylene (PP)
 - Polyvinylidene fluoride (PVDF)


Additional information on the chemical resistance of plastics can be found at www.brand.de.

3 Scope of delivery

Transferpette® S, variable type, marked DE-M, with quality certificate, shelf mount, and silicone grease.

4 Functional and operational components

Front

Adjustable finger rest	<p>The laboratory instrument Transferpette® S has an adjustable finger rest. This allows you to adapt the pipette to your preferred grip, see " Pipetting</p> <p>The instrument can be labeled at the finger rest: to do so, remove the label window from the finger rest and take out the label strip.</p>
Volume-change protection	<p>The volume adjustment lock prevents the volume from being changed during pipetting work. To unlock it, slide the volume adjustment lock toward the pipetting button.</p>
Volume display	<p>The numbers in the display are read from top to bottom; the white dash corresponds to the decimal point.</p>
Counter mechanism protection	<p>Once the volume adjustment lock is released, set the volume using the volume-setting wheel. The integrated counter mechanism protection overrides the volume-setting wheel once the maximum or minimum volume setting is reached: the volume-setting wheel remains rotatable but no longer adjusts the volume.</p>
Hand grip	<p>Screw the pipetting unit into the hand grip. Attach the tip onto the tip cone.</p>

Rear

Permanent adjustment to factory settings:
Easy Calibration

Temporary adjustment to
changing conditions:
User Adjustment

User Adjustment scale

Cover

Seal

Label window

Serial number and
product markings

QR code:

Use your smartphone, tablet,
or webcam to scan and visit
the following website:
www.brand.de/myproduct

The website contains
serialized information about
your laboratory instrument.

Easy Calibration Technology

The Easy Calibration function is located in the finger rest and is used to reset the pipette to its factory settings (see Adjusting – User and factory adjustment, p. 70).

User Adjustment Technology

The hand grip also includes the User Adjustment function. This allows the pipette to be adjusted for specific liquids and dispensing conditions. See Temporary adjustment: User Adjustment, p. 72

The switch for user adjustment is located behind the cover. A seal is applied at delivery. Remove it upon first use and dispose of it.

QR Code and serialized information


The QR code provides access to www.brand.de/myproduct and serialized information about your pipette.

If you wish to access the MyProduct information without the QR code, you will also need the order number and serial number of your pipette.

5 Pipetting

- a. Turn the finger rest to a comfortable working position.

- b. Slide the volume adjustment lock in the direction shown, against slight resistance.


- c. Set the volume using the volume-setting wheel.

NOTICE If the volume-setting wheel is turned beyond the maximum or minimum volume, it glides over the volume adjustment and thus protects the counter mechanism from damage.

- d. Close the volume adjustment lock.
- ⇒ The volume-setting wheel can still be turned but will no longer change the volume.

- e. Attach the pipette tip. Ensure a secure fit.


The 2 ml to 10 ml pipettes should only be used with a built-in PE filter (see UV disinfection, p. 76).

NOTICE Pipette tips are disposable products!

- f. Press the pipetting button down to the first stop.

NOTICE We recommend rinsing the pipette tip five times with the liquid (aspirate and dispense) before pipetting in order to ensure maximum precision and accuracy.

g. Immerse the tip in the liquid.


h. Slowly release the pipetting button.

⇒ Liquid is aspirated.

Keep the tip submerged until the volume is fully aspirated. Increase the wait time when pipetting more viscous liquids or larger volumes.

⚠CAUTION! Do not lay the pipette down with filled tips. Contamination may occur!

Volume range	Immersion depth [mm]	Wait time [s]
0.1–1 μ l	1–2	1
1–100 μ l	2–3	1
100–1000 μ l	2–4	1
> 1,000 μ l	3–6	3

- i. To dispense the liquid, hold the tip at an angle against the vessel wall, slowly press the pipetting button, and wipe the tip.

To improve accuracy, comply with the corresponding wait time for serums, highly-viscous or low-density fluids.

To fully empty the tip, press the pipetting button down to the second stop (Fig. f).

- j. To remove the tip, hold the pipette over a container and press the tip ejection button.

6 Storage

⚠CAUTION! Do not hang the pipette with a filled tip in the holder. Contamination may occur!

You can also hook the Transferpette® S in the holder or rack with an adjustable finger rest.

Rack for 1 pipette

You can use the rack (Item no. [703440](#)) for 1 pipette directly without a holder (stands for different nominal pipette volumes are available).

Rack for 6 pipettes

You can also use the holder for the Transferpette® S in the 6-position rack (Item no. [704807](#)) for the Transferpette® S. To do so, proceed as follows:

- Insert the holder at an angle into the recess and press it downward.
- The holder snaps into place against resistance. Then hang the pipette into the holder.

NOTICE Transferpette® S holders do not secure the Transferpette® S against falling and cannot be used for this pipette.

7 Checking the volume

We recommend testing the device every 3 to 12 months, depending on the level of use. However, the testing cycle can be adapted to meet individual requirements. The complete testing procedure (SOP) can be downloaded at www.brand.de.

You can download the detailed test instructions (SOP) at www.brand.de. For GLP- and ISO-compliant evaluations and documentation, we recommend the EASYCAL™ calibration software from BRAND. A demo version can be downloaded from <https://shop.brand.de/>.

Gravimetric volume testing of the pipette is carried out according to the following steps and complies with DIN EN ISO 8655:2022.

- a. Set User Adjustment to 0 (for procedure, see Temporary adjustment: User Adjustment, p. 72)
- b. Set the maximum specified instrument volume (for procedure, see Pipetting, p. 61).
- c. Condition the pipette before testing by aspirating and dispensing the test liquid (distilled water) with a pipette tip five times.
- d. Aspirate the test liquid and pipette into the weighing vessel.
- e. Weigh the pipetted amount with an analysis scale. (refer to the operating manual of the balance manufacturer.)
- f. Calculate the pipetted volume. In doing so, take into account the temperature of the test liquid.
- g. At least 10 pipetting series and weighings in three volume ranges (100%, 50%, 10%) are recommended. Two tips must be used for each volume range to be tested.

Calculation (for nominal volume)

x_i = weighing results

n = number of weighings

V_0 = nominal volume

Z = Correction factor (e.g. 1.0029 µl/mg at 20°C, 1013 hPa)

Mean:

$$\bar{x} = \frac{\sum x_i}{n}$$

Mean volume:

$$\bar{V} = \bar{x} * Z$$

Accuracy*:

$$A\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Standard deviation*:

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

Coefficient of variation*:

$$CV\% = \frac{100 s}{\bar{V}}$$

*) Accuracy and coefficient of variation are calculated according to the formulas of statistical quality control.

8 Accuracy table

Volume range [μl]	Partial volume [μl]	$A \leq \pm \%$	$CV \leq \pm \%$	Tip type [μl]
0.1-1	1	2	1.2	0.1-20
	0.5	4	2.4	
	0.1	20	12	
0.1-2.5	2.5	1.4	0.7	0.5-20
	1.25	2.5	1.5	
	0.25	12	6	
0.5-10	10	1	0.5	0.5-20
	5	1.6	1	
	1	7	4	
2-20 gray	20	0.8	0.4	0.5-20
	10	1.2	0.7	
	2	5	2	
2-20 yellow	20	0.8	0.4	2-200
	10	1.2	0.7	
	2	5	2	
5-50	50	0.8	0.3	2-200
	25	1.2	0.5	
	5	4	2	
10-100	100	0.6	0.2	2-200
	50	0.8	0.4	
	10	3	1	
20-200	200	0.6	0.2	2-200
	100	0.8	0.3	
	20	3	0.6	
30-300	300	0.6	0.2	5-300
	150	0.8	0.3	
	30	3	0.6	

Volume range [μl]	Partial volume [μl]	A $\leq \pm$ %	CV $\leq \pm$ %	Tip type [μl]
100-1000	1000	0.6	0.2	50-1000
	500	0.8	0.3	
	100	3	0.6	
250-2500	2500	0.6	0.2	500-5000
	1250	0.8	0.3	
	250	3	0.6	
500-5000	5000	0.6	0.2	500-5000
	2500	0.8	0.3	
	500	3	0.6	
1000-10000	10000	0.6	0.2	1000-10000
	5000	0.8	0.3	
	1000	3	0.6	

*A = Accuracy, CV = Coefficient of Variation

Final test values based on the nominal volume (= max. volume) printed on the device and the specified partial volumes at the same temperature (20 °C/68 °F) of the device, surroundings and distilled water, in accordance with DIN EN ISO 8655.

9 Adjusting – User and factory adjustment

You have the following options to adjust the instrument:

- **Factory adjustment:**
The factory adjustment is used for permanent calibration of the instruments to aqueous media according to ISO 8566 in cases of volume deviations.
- **Temporary User Adjustment:**
The User Adjustment is used for temporary volume adjustment under changing conditions. It can be reset to the original state (factory adjustment).

9.1 Permanent factory adjustment: Easy Calibration

The instrument is factory-calibrated to aqueous solutions in accordance with ISO 8655. If it is determined that the pipette is inaccurate, it can be adjusted using the Easy Calibration technique.


a.


- a. Check whether the User Adjustment is set to 0 (see Temporary adjustment: User Adjustment, p. 72).

NOTICE If the User Adjustment is not set to **0**, the pipette will be misaligned when attempting to perform factory adjustment. In this case, set the User Adjustment to **0** and repeat the factory adjustment as described.


- b. Perform a volume check and determine the actual value (see Checking the volume).

- c. Slightly lift and set aside the labeling window (1) on the finger rest with your thumb. Use a paperclip or an unused pipette tip to remove and dispose of the protective foil (2).

- d. Slide the red adjustment slider back completely, lift the volume-setting wheel (decoupling) and release the adjustment slider.

- e. Set the volume adjustment lock to the UNLOCK position and adjust the previously determined actual volume using the volume-setting wheel. Position UNLOCK (see Pipetting, p. 61 > “Set volume”). Set the volume adjustment lock back to the LOCK position. After each adjustment, a volume check is recommended.

- f. Slide the adjustment slider completely back again, allow the volume-setting wheel to slide downward and release the adjustment slider. If the volume-setting wheel does not slide down easily, move it slightly back and forth until it clicks into place. Reinsert the label window.

NOTICE The change to factory settings is indicated by the red adjustment slider now visible in the label window.

9.2 Temporary adjustment: User Adjustment

Temporary User Adjustment improves accuracy under conditions that deviate from the factory settings (aqueous medium, ISO 8655). This enables temporary volume corrections under changing conditions because deviations from water in physical properties, temperature differences between liquid and ambient conditions, specific tip designs, and environmental factors can all affect accuracy.

NOTICE User Adjustment modifies the volume setting by a certain volume offset (e.g. 100 μl : +2 μl = +2%). If the volume setting changes significantly (e.g. 10 μl : +2 μl = +20%), the adjustment value must be recalculated.

Setting the User Adjustment

- Pry off and remove the cover (1) and seal (2) (e.g., using a paperclip). Dispose of the seal.
- Slide the slider (3) down into the recess and hold it there. Use the volume-setting wheel (4) to set the desired User Adjustment value (see below) on the scale. Release the volume-setting wheel and slowly return the slider (3).

NOTICE If the slider is stuck, gently push it back toward the recess (3) and slowly return it again.

- ⇒ The value is set when the user adjustment value aligns with the marking (5).
- Reinsert the cover (1).
 - Verify the adjustment gravimetrically.

Determining the User Adjustment

Example: Pipetting 180 µl with a 20–200 µl pipette

- Perform control weighings on a precision balance and calculate the actual volume:
Actual volume: 178.4 µl
- Calculate the volume to be corrected:
Volume offset: 1.6 µl (= 180 µl – 178.4 µl)

c. Determine and set the User Adjustment value:

For our 200 µl instrument, each dash corresponds to a step value of 0.2 µl (see assignment table). A volume offset of 1.6 µl is added by setting to +8 (= 1.6 µl / 0.2).

$$\text{Actual volume} = \frac{\text{Mean of liquid weights} / (\text{Density of liquid} - \text{Density of air (0.0012 g/ml)})}{\text{Density of liquid} - \text{Density of air (0.0012 g/ml)}}$$

$$\text{Volume offset} = \text{Target volume} - \text{Actual volume}$$

$$\text{User Adjustment value} = \frac{\text{Volume offset}}{\text{Step value}}$$

Assignment table for User Adjustment

The highlighted column [1] indicates the step value for the respective instrument.															
	-25	-20	-15	-10	-5	-1	0	1	5	10	15	20	25	30	35
The step value corresponds to a volume compensation in µl:															
nominal volume [µl]															
1	-0,025	-0,02	-0,015	-0,01	-0,005	-0,001	0	0,001	0,05	0,01	0,015	0,02	0,025	0,03	0,035
2,5	-0,05	-0,04	-0,03	-0,02	-0,01	-0,002	0	0,002	0,01	0,02	0,03	0,04	0,05	0,06	0,07
10	-0,25	-0,2	-0,15	-0,1	-0,05	-0,01	0	0,01	0,05	0,1	0,15	0,2	0,25	0,3	0,35
20	-0,5	-0,4	-0,3	-0,2	-0,1	-0,02	0	0,02	0,1	0,2	0,3	0,4	0,5	0,6	0,7
50	-1,25	-1	-0,75	-0,5	-0,25	-0,05	0	0,05	0,25	0,5	0,75	1	1,25	1,5	1,75
100	-2,5	-2	-1,5	-1	-0,5	-0,1	0	0,1	0,5	1	1,5	2	2,5	3	3,5
200	-5	-4	-3	-2	-1	-0,2	0	0,2	1	2	3	4	5	6	7
300	-6,225	-4,98	-3,735	-2,49	-1,245	-0,249	0	0,249	1,245	2,49	3,735	4,98	6,225	7,47	8,715
1000	-25	-20	-15	-10	-5	-1	0	1	5	10	15	20	25	30	35
1250	-25	-20	-15	-10	-5	-1	0	1	5	10	15	20	25	30	35
2500	-50	-40	-30	-20	-10	-2	0	2	10	20	30	40	50	60	70
5000	-125	-100	-75	-50	-25	-5	0	5	25	50	75	100	125	150	175
10000	-250	-200	-150	-100	-50	-10	0	10	50	100	150	200	250	300	350

← Volume offset for excess volume
Volume offset for missing volume →

NOTICE The table shows the mechanical relationship between the steps of the User Adjustment. The volume changes indicated are approximate values and apply to the entire volume range of the instrument.

User Adjustment calculation tool

www.brand.de/uad

Recognizing user adjustment

If the red switch is visible on the back of the instrument, it has already been adjusted by a user with the user adjustment. Check whether this adjustment still suits your application (e.g., by performing a control weighing of the pipetted volume). Reset User Adjustment if necessary.

Restore factory adjustment, reset User Adjustment

To reset the User Adjustment, set it to 0 on the scale. This restores the factory adjustment state. We recommend performing a volume check afterward.

10 Disinfection/Autoclaving

10.1 Autoclaving

NOTICE Check the effectiveness of autoclaving yourself!

Maximum safety is achieved through vacuum sterilization. We recommend using sterilization bags.

- a. Eject the pipette tip.
- b. Pack the instrument in a sterilization bag, observing any applicable packaging regulations.
- c. Autoclave the complete pipette without further disassembly. Recommendation for autoclaving according to DIN EN 285 (see table below).
- d. Allow the pipette to completely cool and dry.

Temperature:	121 °C
Pressure	2 bar
Holding time in autoclaves	15 min

After autoclaving, tighten the connection between the hand grip and the pipette shaft if necessary.

10.2 UV disinfection

The device is resistant to normal exposure to a UV disinfection lamp. The effects of the UV exposure may cause some color change.

10.3 PE filter

PE filter for Transferpette® S 2 ml, 5 ml + 10 ml:

A hydrophobic PE filter is used to prevent liquid from entering the pipette.

Change the filter if it becomes wet or contaminated.

- a. Use a flat object, such as a screwdriver.
- b. Remove the filter carefully, without damaging the tip cone.

Remove the filter before autoclaving!


The device can also be operated without a filter.

11 Maintenance

- a. Check the pipette tip cone for damage.
- b. Inspect the piston and seal for contamination.
- c. Check the device for leaks.

We recommend using the BRAND leak detector, the BRAND PLT unit. As an alternative to this, aspirate a sample and hold the device vertically for approx. 10 s. If a drop forms at the pipette tip, follow the malfunction remedy, see Troubleshooting, p. 81 .

11.1 Cleaning – volumes up to 1,000 μL


- Detach the pipette shaft (S) from the hand grip by unscrewing it.
- Unscrew the upper part of the ejector unit (A) from the pipette shaft.
- Pull out the shaft (B, C and D) from the lower part of the ejector unit (E).
- Unscrew the piston unit (B).

NOTICE Do not disassemble the piston unit (B) any further!

- Remove the seal with spring (C) (not possible on Transferpette® S 1 μL , 2.5 μL and 10 μL !).
- Clean the parts shown with a soap solution or isopropanol, and then rinse with distilled water.
- Dry the parts (max. 120 °C/248 °F).
- Grease piston and seal with a very thin layer of supplied silicone grease.

Reassemble the cooled parts in reverse order. Only hand-tighten the piston unit and the upper part of the ejector unit (A, B).

11.2 Cleaning – volumes of 2.5, 5, and 10 ml

- Remove the entire shaft (S) from the hand grip by rotating at the upper end of the ejector (F) and remove the filter (K) from the bottom part of the shaft (H).
- Separate the bottom part of the ejector (F') by unscrewing it from the upper part of the ejector (F).
- Unscrew and dismantle the piston unit (G) with the ejector spring (I) and the bottom part of the shaft (H).
- Remove the O-ring-seal from the piston unit and clean it.

NOTICE Do not disassemble the piston unit (G) any further!



- Clean the piston unit (G) and the bottom part of the shaft (H) with a soap solution or isopropanol, and then rinse with distilled water.
- Dry the parts (max. 120 °C/248 °F) and allow them to cool.
- Carefully lubricate the inside and outside of the O-ring (O) and mount it on the piston.

Reassemble the individual components in reverse order.

12 Troubleshooting

Problem	Possible cause	Corrective action
Tip dripping (device leaking)	Unsuitable tip	Only use high-quality tips
	Tip not seated tightly	Firmly press tip on
The instrument does not aspirate or aspirates too little; the dispensed volume is too low	Seal contaminated	Clean seal
	Seal or cone is damaged	Replace seal or shaft
	Piston is contaminated or damaged	Clean or replace piston
Aspiration is very slow	Shaft is clogged	Clean shaft
Pipette is misadjusted	Calibrated with altered User Adjustment	Set User Adjustment to 0. Readjust
Dispensed volume too large	Pipetting button pressed too far (to the over-stroke point) before aspirating	Ensure proper handling.
Piston sluggish	Piston is contaminated or not greased	Clean piston and apply grease

13 Product markings

Symbol or number	Meaning
	Read the user manual.
XXZXXXXX	Serial number
	The instrument is marked in accordance with the German Measurement and Calibration Act as well as the Measurement and Calibration Regulation. DE-M (DE for Germany), framed by a rectangle, as well as the two last digits of the year the marking was affixed.
	Autoclavable up to the temperature shown
Data matrix code or quick response code	The codes link to the BRAND MyProduct website.
www.brand.de/ip	Hyperlink to the BRAND patent page

14 Order Information

Various pipettes

Nominal volume	Cat. No.
0.1–1 µl	706868
0.1 - 2.5 µl	706869
0.5 - 10 µl	706870
2 - 20 µl (gray)	706871
2 - 20 µl (yellow)	706872
5 - 50 µl	706873
10 - 100 µl	706874
20 - 200 µl	706878
30 - 300 µl	706879
100 - 1,000 µl	706880
250 - 2500 µl	706881
500 - 5,000 µl	706882
1,000 - 10,000 µl	706884

Accessories

Table stand for
1 pipette
Cat. No. [703440](#) or
[705384](#)

Wall mount
Cat. No. [704882](#)

Table stand for 6 single-
channel or multi-
channel pipettes
(Can be used with the
holders of the
Transferpette® pro)
Cat. No. [704807](#)

Shelf mount
Cat. No. [704881](#)

14 Order Information

Labeling window
Cat. No. [704752](#)

Labeling foil
Cat. No. [704753](#)

Filter for volume range 2
- 5 ml, PU 25 pcs
Cat. No. [704652](#)

Silicone grease
volume range up to
1000 μ l
Cat. No. [705502](#)


Filter for volume range
up to 10 ml, PU 25 pcs
Cat. No. [704653](#)

PLT unit (pipette leak
detector)
Cat. No. [703970](#)


Spare parts - volumes up to 1000 μl 

Vol.	A	B	C	D	E
0.1– 1 μl	705513	704600	—	704718 *	704705
0.1– 2.5 μl	705513	704667	—	704717	704706
0.5– 10 μl	705513	704601	—	704721 *	704707
2 – 20 μl (gray)	705513	704602	704610	704727	704710
2 – 20 μl (yellow)	705513	704602	704610	704723	704710
5– 50 μl	705513	704615	704617	704722	704711
10– 100 μl	705513	704654	704661	704724	704712
20– 200 μl	705513	704655	704662	704725	704713
30 – 300 μl	705513	704668	704664	704729	704714
100– 1000 μl	705513	704656	704663	704726	704715

* Seal permanently installed in shaft – not removable!

NOTICE The appearance and dimensions of the spare parts correspond to the respective nominal volume.

Spares - volumes 2.5, 5 and 10 ml

Vol.	F + F'	G	H	I	O
0.25– 2.5 ml	704755	704669	704689	704626	7228
0.5– 5 ml	704756	704606	703247	704626	7228
1– 10 ml	704757	704607	704628	704626	7228

NOTICE The appearance and dimensions of the spare parts correspond to the respective nominal volume.

15 Repairs

15.1 Sending for repair

NOTICE

The transportation of hazardous materials without permission is prohibited by law.

Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

Within USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

Contact addresses

Germany:
BRAND GMBH + CO KG

USA and Canada:
BrandTech® Scientific, Inc.

15 Repairs

Otto-Schott-Straße 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T +91 22 42957790
F +91 22 42957791
info@brand.co.in
www.brand.co.in

11 Bokum Road
Essex, CT 06426-1506 (USA)
T +1-860-767 2562
F +1-860-767 2563
info@brandtech.com
www.brandtech.com

China:

BRAND (Shanghai) Trading Co., Ltd.
Rm 201-202, North Tower,
No. 199 Kaibin Rd, Xuhui District, Shanghai
Shanghai 200030 (P.R. China)
T +86 21 6422 2318
F +86 21 6422 2268
info@brand.com.cn
www.brand.cn.com

16 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring devices. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the instrument. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive liquids.

The detailed testing instructions are available for download on www.brand.de or www.brandtech.com.

BRAND also offers the possibility to have your instruments calibrated by our factory calibration service or our accredited calibration laboratory. Simply send us the instruments to be calibrated, accompanied by details about which type of calibration you would like. The instruments will be returned to you after a few days. The instruments are supplied with a detailed calibration certificate or a calibration slip in accordance with DIN EN ISO/IEC 17025. More information can be obtained from your dealer or directly from BRAND. The order document is available for download on www.brand.de (see Service and Support).

For customers outside Germany

If you would like to use our calibration service, please contact one of our service partners in your region. If a factory calibration is requested, the service partner can forward the instrument to BRAND.

17 Information about your laboratory device

The online service MyProduct (<https://www.brand.de/myproduct>) offers quality certificates, accessories, and technical documentation for your laboratory instrument Transferpette® S. By entering the serial and article number there, you will receive information specific to your individual instrument.

On Transferpette® S, you will also find serialized information encoded in the quick response code. This also links to the MyProduct website, where you can access this operating manual and additional certificates for your instrument.

You will also find a data matrix code on some instruments (Transferpette® S, HandyStep® touch, and HandyStep touch® S). Scan this using a standard reader app to access the above-mentioned information via the URL <https://www.brand.de/myproduct>.

18 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

USA and Canada:

Find more warranty information on www.brandtech.com.

19 Disposal

Before disposal, observe the relevant national disposal regulations and ensure that the product is disposed of properly.

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